



Rotational Viscosity Standards

Product Code	Type	mPa.s (cP) 20°C	mPa.s (cP) 25°C	Pack Size
VIS-RT5-600	RT5	5.431	4.972	600mL
VIS-RT10-600	RT10	10.89	9.972	600mL
VIS-RT50-600	RT50	54.43	49.53	600mL
VIS-RT75-600	RT75	84.55	77.03	600mL
VIS-RT100-600	RT100	111.6	101.5	600mL
VIS-RT250-600	RT250	275.7	250.7	600mL
VIS-RT350-600	RT350	386.6	351.5	600mL
VIS-RT500-600	RT500	561.1	509.4	600mL
VIS-RT1K-600	RT1000	1126	1018	600mL
VIS-RT5K-600	RT5000	5492	4962	600mL
VIS-RT10K-600	RT10000	11142	10040	600mL
VIS-RT12K-600	RT12000	14159	12769	600mL
VIS-RT30K-600	RT30000	33168	29907	600mL
VIS-RT60K-600	RT60000	67226	60586	600mL
VIS-RT100K-600	RT100000	110159	99565	600mL

All values stated are nominal - 08/12/2008