

TECHLAB

Matériaux de Référence Certifiés & Equipements de Laboratoire

Plastic Additives



TECHLAB

Introduction

Plastics and other polymeric materials have become indispensable in our everyday lives. Although they offer many benefits, hazardous chemicals may be present in these materials. These hazardous materials can be introduced either intentionally as additives, or unintentionally as pollutants.

We have collected or synthesized many of these polymer adjuncts and are offering them as Certified Reference Materials in this guide/catalog. Calibrating with certified standards adds an additional layer of confidence in the analysis that can aid in complying with government/industry regulations.

The occurrence, toxicity and analytical methods used in the detection and monitoring (for both presence and levels) of these chemicals are more thoroughly described in the book "Handbook for the Chemical Analysis of Plastic and Polymer Additives" 2nd Edition by CRC Press. Both manufacturers and distributors of plastic and polymeric materials will find the CRC book to be an authoritative source of information that complements this guide.

Both this guide and the book organize additives by the chemical function, e.g. Antioxidant. Manufacturers can easily find standards that match their particular application for the following product categories:

Medical devices	Pharmaceutical packaging	Wire and cable
Food packaging	Toys	

Below is a listing of regulations that require analysis for many of these additives.

- EU Directive 2002/96/EC WEEE (Waste Electrical and Electronic Equipment) establishes limits for the content of a product that must be recyclable or reusable.
- EU Directive 2003/11/EC ROHS (Restriction of the use of certain hazardous substances) restricting the use of six toxins from most electronic and electrical equipment.
- EU Directive 2002/72/EC relating to plastic materials and articles intended to come in contact with foodstuffs.
- EU Directive 2015/863 restricting the use of four phthalates.
- Annex XVII of REACH, annex 51, restriction of phthalates in children's toys and child care products.
- Rule 82 FR 49938 Consumer Product Safety Commission (CPSC) ban of five phthalates in children's toys and child care products.
- United States Environmental Protection Agency (USEPA) – Methods 506, 606 and 8061 regulating phthalates and adipates.

Phthalates

Recently, products containing phthalates have received intense scrutiny. Phthalates are used primarily as plasticizers. Plasticizers lower the glass transition temperature of a plastic/polymer, and impart flexibility, durability and longevity to these types of materials by acting as softening agents.

Due to their low cost, versatility and effectiveness, they are end-use components in products such as electronics, paints, adhesives, building materials and toys. Lower molecular weight phthalates are easily released into the environment because there is no chemical bond between the phthalate and the plastic/polymer matrix. Leaching and atmospheric release of these compounds increases as the substrate ages and/or weathers, and has resulted in phthalates becoming a major environmental contaminant. This is important because phthalates are considered to be potential endocrine-disrupting agents with human exposure through direct contact, ingestion or inhalation.

Concern over the adverse health effects has prompted regulatory changes and led to a permanent ban of these plasticizers in child-care products and toys (see above).

This prohibition has introduced the use of phthalate replacements.

Plastic Additives

Table of Contents	
Phthalates	2-3
Isophthalates	4
Terephthalates	4
Monophthalates	4-5
Deuterated Phthalates	6
Phthalates - Industrial	7
Phthalate Replacements	7-8
Phthalate EPA Methods	9
Accelerants	10
Antifoams	11
Antidegradants	11
Antioxidants	12-19
Antiozonates	20
Blowing Agents	21
Coupling Agents	21
Cross Linking Agents	22
Flame Retardants	23-25
Plasticizers	26-33
Processing Aids	34
Retarders	34-35
Stearates	35
UV Stabilizers	36-37
Vegetable Oils	37
Bisphenol Analogs	38-39
Plastic Packaging Testing	40
Dyes and Breakdown Products	41
Index	

**Custom Quotation
Request Form**

Company :

Contact : Job Function:

Address :

Postal Code: Town: Country:

Telephone : FAX :

Email : URL :

Area of Interest :

Product Description :

Concentration :

Matrix :

Concentration Units :

- ng/ml
- µg/ml
- mg/ml
- wt. %
- vol. %

Requested Quantity :**Organic** 5 x 1 ml**Inorganic** ____ x 500 ml Others ____ x ____**Component(s)**

1.
2.
3.
4.
5.
6.
7.
8.
9.
10.
11.
12.
13.
14.
15.
16.

Concentration (if various)


-
-
-
-
-
-
-
-
-
-
-
-
-
-
-
-

Phthalates Table of Contents

Phthalates	2-3
Isophthalates	4
Terephthalates	4
Monophthalates	4-5
Deuterated Phthalates	6
Phthalates - Industrial	7
Phthalate Replacements	8
Phthalate EPA Methods	9

The perfect companion for your analysis!

This reference book contains the compounds in this catalog, with important reference data to aid in testing and compliance.



Each Compound has:

Chemical Information

- Structure
- CAS Number (where applicable)
- RTECS Number (where available)
- Formula
- Molecular Weight
- IUPAC Name, other common names and some popular brand names

Other Important Information

- Application
- Regulatory
- Environmental Impact
- Point of Release
- Toxicological Data

Analytical Data

- Mass Spectrum with key ions tabulated
- Chromatogram with conditions

Physical Properties

- Appearance
- Melting and Boiling Points
- Stability
- Solubilities in several common solvents

Cat. No: BOOK-PLAS-002

Handbook for the Chemical Analysis of Plastic and Polymer Additives, 2nd Ed.

The Second Edition of this handbook provides the necessary tools for chemists to obtain a more complete listing of additives present in a particular polymeric matrix. This edition features:


- Updated material to include the most recent additives available
- Actual analytical data for each chemical along with the description and methods used for obtaining the results
- Toxicological and environmental impact of each product highlighted
- Summaries of regulatory and health information in a convenient "one-step" format

With 50 additional compounds, this 2nd edition nearly doubles the number of additives in several categories including processing aids, anti-static compounds, mold release products and blowing agents. It includes a listing that can be cross-referenced by trade name, chemical name, CAS number and even key mass unit ions from the GC/MS run. Also included are case studies related to "real-world" issues, tips for analysis in challenging matrices and more.

Phthalates

Phthalates


Benzyl butyl phthalate



**CAS 85-68-7 MF C₁₉H₂₀O₄ MW 312.26 PS L
SG 1.13 g/cm³ MP N/A BP 370-380 °C
FP 198 °C**

Matrix	Cat. No.	Unit
NEAT	ALR-082N	100 mg
100 µg/mL in MeOH	ALR-082S	1 mL
5 mg/mL in MeOH	AS-E0065	1 mL


bis(2-Methoxyethyl)phthalate



**CAS 117-82-8 MF C₁₄H₁₈O₆ MW 282.29 PS L
SG 1.17 g/cm³ MP N/A BP 340 °C FP 185 °C**

Matrix	Cat. No.	Unit
NEAT	J-106	100 mg


Diamyl phthalate



**CAS 131-18-0 MF C₁₈H₂₆O₄ MW 306.40 PS L
SG 1.03 g/cm³ BP 342 °C FP 118 °C**

Matrix	Cat. No.	Unit
NEAT	ALR-098N	100 mg
100 µg/mL in MeOH	ALR-098S	1 mL


bis(4-Methylpentyl)phthalate



**CAS 71850-09-4 MF C₂₀H₃₀O₄ MW 334.45 PS S
SG N/A BP N/A FP N/A**

Matrix	Cat. No.	Unit
NEAT	PHTH-022N	100 mg
100 µg/mL in MeOH	PHTH-022S	1 mL


Dibenzyl phthalate



**CAS 523-31-9 MF C₂₂H₁₈O₄ MW 346.38 PS S
SG 1.25 g/cm³ MP 40-42 °C BP >400 °C
FP >150 °C**

Matrix	Cat. No.	Unit
NEAT	J-104	100 mg


bis(2-n-Butoxyethyl)phthalate



**CAS 117-83-9 MF C₂₀H₃₀O₆ MW 366.45 PS L
SG 1.06 g/cm³ MP N/A BP 270 °C FP 205 °C**

Matrix	Cat. No.	Unit
NEAT	J-112	100 mg


bis(2-Ethoxyethyl)phthalate



**CAS 605-54-9 MF C₁₆H₂₂O₆ MW 310.34 PS S
SG 1.12 g/cm³ MP 34 °C BP 345 °C FP N/A**

Matrix	Cat. No.	Unit
NEAT	J-111	100 mg


bis(2-Ethylhexyl)phthalate (DEHP)



**CAS 117-81-7 MF C₂₄H₃₈O₄ MW 390.56 PS L
SG 0.98 g/cm³ MP N/A BP 385 °C FP 218 °C**

Matrix	Cat. No.	Unit
NEAT	ALR-097N	100 mg
100 µg/mL in MeOH	ALR-097S	1 mL
1 mg/mL in MeOH	APP-9-029-10X	1 mL


2-Butoxy-2-oxoethyl butyl phthalate



**CAS 85-70-1 MF C₁₈H₂₄O₆ MW 336.38 PS L
SG 1.10 g/cm³ MP N/A BP 345 °C FP 199 °C**

Matrix	Cat. No.	Unit
NEAT	J-115	100 mg


Diallyl phthalate



**CAS 131-17-9 MF C₁₄H₁₄O₄ MW 246.26 PS L
SG 1.11 g/cm³ MP N/A BP 165 °C FP 165 °C**

Matrix	Cat. No.	Unit
NEAT	J-002	100 mg

Dicyclohexyl phthalate




**CAS 84-61-7 MF C₂₀H₂₆O₄ MW 330.42 PS S
SG 1.14 g/cm³ MP 61-66 °C BP 235 °C
FP 207 °C**

Matrix	Cat. No.	Unit
NEAT	J-004	100 mg
100 µg/mL in MeOH	ALR-099S	1 mL
1 mg/mL in AcCN	AS-E0318	1 mL

Phthalates

Phthalates


Didodecyl phthalate



CAS 2432-90-8 MF C₃₂H₅₄O₄ MW 502.77
PS L or S SG 1.05 g/cm³ MP 21-23 °C BP 256 °C
FP >200 °C

Matrix	Cat. No.	Unit
NEAT	PHTH-018N	100 mg
100 µg/mL in MeOH	PHTH-018S	1 mL


Diisobutyl phthalate



CAS 84-69-5 MF C₁₆H₂₂O₄ MW 278.34 PS L
SG 1.04 g/cm³ MP N/A BP 327 °C FP 180 °C

Matrix	Cat. No.	Unit
NEAT	J-113	100 mg


Di-n-octyl phthalate



CAS 117-84-0 MF C₂₄H₃₈O₄ MW 390.56 PS L
SG 0.98 g/cm³ MP N/A BP 109 °C

Matrix	Cat. No.	Unit
NEAT	J-011	100 mg
100 µg/mL in MeOH	ALR-105S	1 mL
5 mg/mL in MeOH	AS-E0067	1 mL


Diethyl phthalate



CAS 84-66-2 MF C₁₂H₁₄O₄ MW 222.24 PS L
SG 1.12 g/cm³ MP N/A BP 298-299 °C FP 160 °C

Matrix	Cat. No.	Unit
NEAT	J-005	100 mg
100 µg/mL in MeOH	APP-9-081	1 mL
1 mg/mL in MeOH	APP-9-081-10X	1 mL
5 mg/mL in MeOH	AS-E0068	1 mL


Diisopentyl phthalate



CAS 605-50-5 MF C₁₈H₂₆O₄ MW 306.40 PS L
SG 1.03 g/cm³ FP 167 °C

Matrix	Cat. No.	Unit
NEAT	J-127	100 mg


Diphenyl phthalate



CAS 84-62-8 MF C₂₀H₁₄O₄ MW 318.32 PS S
SG 1.24 g/cm³ MP 74-76 °C BP 255 °C FP 256 °C

Matrix	Cat. No.	Unit
NEAT	J-013	100 mg


Di-n-heptyl phthalate



CAS 3648-21-3 MF C₂₂H₃₄O₄ MW 362.50 PS L
SG 0.99 g/cm³ MP N/A BP 360 °C FP 113 °C

Matrix	Cat. No.	Unit
NEAT	PHTH-020N	100 mg
100 µg/mL in MeOH	PHTH-020S	1 mL


Diisopropyl phthalate



CAS 605-45-8 MF C₁₄H₁₈O₄ MW 250.29 PS L
MP N/A

Matrix	Cat. No.	Unit
NEAT	PHTH-019N	100 mg


Dihexyl phthalate



CAS 84-75-3 MF C₂₀H₃₀O₄ MW 334.45 PS L
SG 1.01 g/cm³ BP 333 °C FP 177 °C

Matrix	Cat. No.	Unit
NEAT	ALR-100N	100 mg
100 µg/mL in MeOH	ALR-100S	1 mL

Dimethyl phthalate



CAS 131-11-3 MF C₁₀H₁₀O₄ MW 194.18 PS L
SG 1.19 g/cm³ MP N/A BP 282-284 °C
FP 156 °C

Matrix	Cat. No.	Unit
NEAT	J-010	100 mg
100 µg/mL in MeOH	APP-9-088	1 mL
1 mg/mL in MeOH	APP-9-088-10X	1 mL
5 mg/mL in MeOH	AS-E0069	1 mL
0.1 mg/mL in EtOAc	M-8032-IS	1 mL


Property Key

CAS	Chemical Abstract Service Number	SG	Specific Gravity (g/cm ³)
MF	Molecular Formula	MP	Melting Point (°C)
MW	Molecular Weight	BP	Boiling Point (°C)
PS	Physical State (Solid, Liquid)	FP	Flash Point (°C)

Iso, Tere, and Mono Phthalates

Isophthalates


Dimethyl isophthalate



CAS 1459-93-4 MF C₁₀H₁₀O₄ MW 194.18 PS S
SG 1.18 g/cm³ MP 64-68 °C BP 282-285 °C
FP 138 °C

Matrix	Cat. No.	Unit
NEAT	J-009	100 mg

Diphenyl isophthalate




CAS 744-45-6 MF C₂₀H₁₄O₄ MW 318.32 PS S
SG 1.24 g/cm³ MP 136-138 °C FP 256 °C

Matrix	Cat. No.	Unit
NEAT	J-012	100 mg

Terephthalates


bis(2-Ethylhexyl) terephthalate



CAS 6422-86-2 MF C₂₄H₃₈O₄ MW 390.56 PS L
SG 0.99 g/cm³ MP N/A BP 400 °C FP 238 °C

Matrix	Cat. No.	Unit
NEAT	J-121	100 mg


Diethyl terephthalate



CAS 636-09-9 MF C₁₂H₁₄O₄ MW 222.24 PS S
SG 1.15 g/cm³ MP 43-47 °C BP 142 °C
FP >150 °C

Matrix	Cat. No.	Unit
NEAT	J-123	100 mg

Dimethyl terephthalate




CAS 120-61-6 MF C₁₀H₁₀O₄ MW 194.18 PS S
SG 1.36 g/cm³ MP 139-141 °C BP 288 °C
FP 151 °C

Matrix	Cat. No.	Unit
NEAT	J-101	100 mg

Monophthalates

Mono-phthalate esters are the primary phthalate metabolites formed via hydrolysis of one ester bond. It is these compounds that are thought to be toxic agents and are receiving interest as a possible human health issue. Studies have shown that they can produce estrogenic and immune-suppressive effects in humans.


Monobenzyl phthalate (mBzP)



CAS 2528-16-7 MF C₁₅H₁₂O₄ MW 256.25 PS S
SG 1.28 g/cm³ MP 93-94 °C BP 441 °C FP 168 °C

Matrix	Cat. No.	Unit
NEAT	ALR-134N	100 mg
100 µg/mL in AcCN	ALR-134S-CN	1 mL


Monobutyl phthalate (mBP)



CAS 131-70-4 MF C₁₂H₁₄O₄ MW 222.24 PS S
SG 1.17 g/cm³ MP 65-66 °C BP 350-354 °C
FP 138 °C

Matrix	Cat. No.	Unit
NEAT	ALR-135N	100 mg
100 µg/mL in AcCN	ALR-135S-CN	1 mL

Monocyclohexyl phthalate (mBP)



CAS 7517-36-4 MF C₁₄H₁₆O₄ MW 248.27 PS L
SG 1.24 g/cm³ MP N/A BP 410 °C FP 154 °C

Matrix	Cat. No.	Unit
NEAT	ALR-178N	100 mg
100 µg/mL in AcCN	ALR-178S-CN	1 mL

Technical Note


13 mono-phthalates including the monoethylhexyl (mEHP) which is the metabolite of the plasticizer with the greatest yearly production and use on a global basis.

Monophthalates continued on next page

Monophthalates

Monophthalates (continued)


Monoethyl phthalate (mEP)



CAS 2306-33-4 MF C₁₀H₁₀O₄ MW 194.18 PS S SG 1.24 g/cm³ MP 46-48 °C BP 339 °C FP 135 °C

Matrix	Cat. No.	Unit
NEAT	ALR-137N	100 mg
100 µg/mL in AcCN	ALR-137S-CN	1 mL


Monoisobutyl phthalate



CAS 30833-53-5 MF C₁₂H₁₄O₄ MW 222.24 PS L SG 1.17 g/cm³ MP N/A BP 356-357 °C FP 135 °C

Matrix	Cat. No.	Unit
NEAT	ALR-176N	100 mg
100 µg/mL in AcCN	ALR-176S-CN	1 mL


Monooctyl phthalate



CAS 5393-19-1 MF C₁₆H₂₂O₄ MW 278.34 PS L SG 1.09 g/cm³ MP N/A BP 400 °C FP 146 °C

Matrix	Cat. No.	Unit
NEAT	ALR-141N	100 mg
100 µg/mL in AcCN	ALR-141S-CN	1 mL


Mono-2-ethylhexyl phthalate (mEHP)



CAS 4376-20-9 MF C₁₆H₂₂O₄ MW 278.34 PS L SG 1.09 g/cm³ MP N/A BP 408 °C FP 144 °C

Matrix	Cat. No.	Unit
NEAT	ALR-138N	100 mg
100 µg/mL in AcCN	ALR-138S-CN	1 mL


Monoisononyl phthalate



CAS N/A MF C₁₇H₂₄O₄ MW 292.37 PS S MP 37-38 °C

Matrix	Cat. No.	Unit
NEAT	ALR-142N	100 mg
100 µg/mL in AcCN	ALR-142S-CN	1 mL


Mono-n-pentyl phthalate



CAS 24539-56-8 MF C₁₃H₁₆O₄ MW 236.26 PS L SG 1.48 g/cm³ MP N/A BP 376-377 °C FP 140 °C

Matrix	Cat. No.	Unit
NEAT	ALR-177N	100 mg
100 µg/mL in AcCN	ALR-177S-CN	1 mL


Mono-2-heptyl phthalate



CAS N/A MF C₁₅H₂₀O₄ MW 264.32 PS L MP N/A

Matrix	Cat. No.	Unit
NEAT	ALR-143N	100 mg
100 µg/mL in AcCN	ALR-143S-CN	1 mL


Monoisopropyl phthalate



CAS 35118-50-4 MF C₁₁H₁₂O₄ MW 208.21 PS S SG 1.20 g/cm³ MP 100-104 °C BP 343-344 °C FP 133 °C

Matrix	Cat. No.	Unit
NEAT	ALR-179N	100 mg
100 µg/mL in AcCN	ALR-179S-CN	1 mL


Monohexyl phthalate



CAS 24539-57-9 MF C₁₄H₁₈O₄ MW 250.29 PS L SG 1.12 g/cm³ MP N/A BP 375-380 °C FP 142 °C

Matrix	Cat. No.	Unit
NEAT	ALR-175N	100 mg
100 µg/mL in AcCN	ALR-175S-CN	1 mL

Monomethyl phthalate



CAS 4376-18-5 MF C₉H₈O₄ MW 180.16 PS S SG 1.29 g/cm³ MP 81-84 °C BP 328 °C FP 135 °C


Matrix	Cat. No.	Unit
NEAT	ALR-139N	100 mg
100 µg/mL in AcCN	ALR-139S-CN	1 mL

Deuterated Phthalates

Deuterated Phthalates

Eleven deuterated phthalates which can be used as internal standards. To simplify the ordering process, the native and the corresponding deuterated compound are packaged as sets at a reduced price


Dibenzylphthalate-3,4,5,6-d₄



CAS 1015854-62-2 MF C₂₂H₁₄D₄O₄ MW 350.40
PS S MP 40-42 °C BP 276-278 °C

Matrix	Cat. No.	Unit
NEAT only as Set	PHTH-D4-001N	5 mg
100 µg/mL in MeOH	PHTH-D4-001S	1 mL


Dimethyl phthalate-3,4,5,6-d₄



CAS 93951-89-4 MF C₁₀H₆D₄O₄ MW 198.21 PS L SG 1.20 g/cm³ MP N/A BP 282 °C FP 295 °C

Matrix	Cat. No.	Unit
NEAT only as Set	PHTH-D4-007N	5 mg
100 µg/mL in MeOH	PHTH-D4-007S	1 mL


Di-n-pentyl phthalate-3,4,5,6-d₄



CAS 358730-89-9 MF C₁₈H₂₂D₄O₄ MW 310.42 PS L MP N/A BP 342 °C FP 118 °C

Matrix	Cat. No.	Unit
NEAT only as Set	PHTH-D4-009N	5 mg
100 µg/mL in MeOH	PHTH-D4-009S	1 mL


Dicyclohexyl phthalate-3,4,5,6-d₄



CAS 358731-25-6 MF C₂₀H₂₂D₄O₄ MW 334.44 PS S SG 1.16 g/cm³ MP 65-67 °C FP 207 °C BP 426 °C

Matrix	Cat. No.	Unit
NEAT only as Set	PHTH-D4-004N	5 mg
100 µg/mL in MeOH	PHTH-D4-004S	1 mL


Di-n-butyl phthalate-3,4,5,6-d₄



CAS 93952-11-5 MF D₁₆H₁₈D₄O₄ MW 282.37 PS L SG 1.07 g/cm³ MP N/A BP 336 °C FP 340 °C

Matrix	Cat. No.	Unit
NEAT only as Set	PHTH-D4-002N	5 mg
100 µg/mL in MeOH	PHTH-D4-002S	1 mL


Di-n-propyl phthalate-3,4,5,6-d₄



CAS 358731-29-0 MF C₁₄H₁₄D₄O₄ MW 254.31 PS L SG 1.08 g/cm³ MP N/A BP 317-318 °C FP >110 °C

Matrix	Cat. No.	Unit
NEAT only as Set	PHTH-D4-010N	5 mg
100 µg/mL in MeOH	PHTH-D4-010S	1 mL


Diethyl phthalate-3,4,5,6-d₄



CAS 93952-12-6 MF C₁₂H₁₀D₄O₄ MW 226.26 PS L SG 1.14 g/cm³ MP N/A FP 160 °C BP 298-299 °C

Matrix	Cat. No.	Unit
NEAT only as Set	PHTH-D4-005N	5 mg
100 µg/mL in MeOH	PHTH-D4-005S	1 mL


Diisobutyl phthalate-3,4,5,6-d₄



CAS 358730-88-8 MF C₁₆H₁₈D₄O₄ MW 282.37 PS L MP N/A BP 327 °C FP 109 °C

Matrix	Cat. No.	Unit
NEAT only as Set	PHTH-D4-003N	5 mg
100 µg/mL in MeOH	PHTH-D4-003S	1 mL


Di-n-octyl phthalate-3,4,5,6-d₄



CAS 93952-13-7 MF C₂₄H₃₄D₄O₄ MW 394.58 PS L SG 0.96 g/cm³ MP N/A BP 384 °C FP 405 °C

Matrix	Cat. No.	Unit
NEAT only as Set	PHTH-D4-008N	5 mg
100 µg/mL in MeOH	PHTH-D4-008S	1 mL

bis(2-Ethylhexyl) phthalate-3,4,5,6-d₄



CAS 93951-87-2 MF C₂₄H₃₄D₄O₄ MW 394.58 PS L SG 0.98 g/cm³ MP N/A BP 384 °C FP 207 °C

Matrix	Cat. No.	Unit
NEAT only as Set	PHTH-D4-011N	5 mg
100 µg/mL in MeOH	PHTH-D4-011S	1 mL

Set includes 11 Deuterated Phthalates

NEAT Set	PHTH-D4N-SET	11 x 5 mg
SOLUTION Set	PHTH-D4S-SET	11 x 1 mL

Other compounds are available.
contact our Technical Service if you
require additional deuterated or other
labeled compounds.

Industrial and Phthalates Replacements

Industrial Phthalates

The high molecular weight (HMW) phthalates have more than six carbons in the backbone and are synthesized from phthalic acid and mixtures of C9 and C10 alcohols. The two major HMW products are diisononyl phthalate (DINP) and diisodecyl phthalate (DIDP).

Attention has now turned to the analysis of these compounds as they are becoming major players in the plasticizer marketplace. However, due to the synthesis process, GC separation of DINP and DIDP results in a cluster of peaks corresponding to different isomers. Consequently, different analytical approaches based on soft ionization techniques and MS detection have been documented in the literature (1). These new approaches can provide another tool to scrutinize the amounts, environmental fate and potential health effects of these HMW plasticizers.

1. David, F., Sandra, P. and Hancock, P., *Current Trends in Mass Spectrometry*, May 2011)

Compound	CAS No.	Conc.	Matrix	Cat. No.	Unit
Benzyl 2-ethylhexyl phthalate	27215-22-1	100 mg	NEAT	ALR-165N	100 mg
		100 µg/mL	MeOH	ALR-165S	1 mL
n-Butyl benzyl phthalate	85-68-7	10 mg	NEAT	PHTH-014N	10 mg
		100 µg/mL	MeOH	PHTH-014S	1 mL
Butyl cyclohexyl phthalate (Tech mix)	84-64-0	100 mg	NEAT	J-122	100 mg
n-Butyl isobutyl phthalate (Tech mix)		10 mg	NEAT	PHTH-013N	10 mg
		100 µg/mL	MeOH	PHTH-013S	1 mL
Butyl octyl phthalate (Tech mix)	84-78-6	100 mg	NEAT	J-001	100 mg
Decyl octyl phthalate (Tech mix)		10 mg	NEAT	PHTH-012N	10 mg
		100 µg/mL	MeOH	PHTH-012S	1 mL
Didecyl phthalate	84-77-5	100 mg	NEAT	J-120	100 mg
Diisodecyl phthalate (Tech mix)	26761-40-0	100 mg	NEAT	ALR-101N	100 mg
		100 µg/mL	MeOH	ALR-101S	1 mL
Diisoheptyl phthalate	71888-89-6	100 mg	NEAT	PHTH-017N	100 mg
		100 µg/mL	MeOH	PHTH-017S	1 mL
Diisohexyl phthalate (Tech mix)	68515-50-4	100 mg	NEAT	J-007	100 mg
Diisononyl phthalate (Tech mix) [C8 to C10 Isomers]	68515-48-0	100 mg	NEAT	ALR-102N	100 mg
		100 µg/mL	MeOH	ALR-102S	1 mL
Diisooctyl phthalate (Tech mix) [C8 Isomers]	27554-26-3	100 mg	NEAT	ALR-103N	100 mg
		100 µg/mL	MeOH	ALR-103S	1 mL
Dinonyl phthalate	84-76-4	100 mg	NEAT	J-105	100 mg
Hexyl 2-ethylhexyl phthalate (Tech mix)	75673-16-4	100 mg	NEAT	J-016	100 mg
Isobutyl benzyl phthalate (Tech mix)		10 mg	NEAT	PHTH-015N	10 mg
		100 µg/mL	MeOH	PHTH-015S	1 mL
Isobutyl cyclohexyl phthalate (Tech mix)	5334-09-8	100 mg	NEAT	J-014	100 mg
n-Octyl n-decyl phthalate (Tech mix)	119-07-3	100 mg	NEAT	J-015	100 mg
Pentyl isopentyl phthalate (Tech mix)	84777-06-0	10 mg	NEAT	PHTH-016N	10 mg
		100 µg/mL	MeOH	PHTH-016S	1 mL

Phthalate Replacements

World-wide concern over environmental and health-related factors associated with phthalates has led to restrictions of use in a wide array of products. This has resulted in the plastics industry generating a variety of alternatives.

In response, we have developed a phthalate replacement product line comprised of 42 compounds representing 18 chemical classes.

Compound	CAS No.	Concentration	Cat. No.	Unit
Azelaic Acid Derivatives				
Diisodecyl azelate	28472-97-1	1000 µg/mL in Acetone	PLAS-PL-075S-A	1 mL
Diisooctyl azelate	26544-17-2	1000 µg/mL in Acetone	PLAS-PL-076S-A	1 mL
Dimethyl azelate	1732-10-1	1000 µg/mL in Acetone	PLAS-PL-077S-A	1 mL
Di-n-hexyl azelate	109-31-9	1000 µg/mL in Acetone	PLAS-PL-078S-A	1 mL
Di(2-ethyl hexyl) azelate	103-24-2	1000 µg/mL in Acetone	PLAS-PL-081S-A	1 mL
Adipic Acid Derivatives				
Di(tridecyl) adipate	16958-92-2	1000 µg/mL in Acetone	PLAS-PL-079S-A	1 mL
Di(n-heptyl, n-nonyl) adipate	68515-75-3	1000 µg/mL in Hexane	PLAS-PL-080S	1 mL
Diisobutyl adipate	141-04-8	1000 µg/mL in Hexane	PLAS-PL-082S	1 mL
Diisodecyl adipate	27178-16-1	1000 µg/mL in Hexane	PLAS-PL-083S	1 mL
Dimer Acid Derivatives				
bis(2-Hydroxyethyl) dimerate	68855-78-7	1000 µg/mL in Hexane	PLAS-PL-084S	1 mL
Epoxy Derivatives				
Epoxidized linseed oil	8016-11-3	1000 µg/mL in Toluene	PLAS-PL-085S-T	1 mL
2-Ethylhexyl epoxy tallate	61789-01-3	1000 µg/mL in Hexane	PLAS-PL-086S	1 mL
Fumaric Acid Derivative				
Dibutyl fumarate	105-75-9	1000 µg/mL in Hexane	PLAS-PL-087S	1 mL
Glycerol Derivative				
Glycerol triacetate	102-76-1	1000 µg/mL in Hexane	PLAS-PL-088S	1 mL
Isobutyrate Derivative				
2,2,4-Trimethyl-1,3-pentanediol-diisobutyrate	6846-50-0	1000 µg/mL in Hexane	PLAS-PL-089S	1 mL

See Plasticizer section,
pages 26-33, for
structures, molecular
formulas and molecular
weights.

Replacement Phthalates continued on next page

Phthalate Replacements (cont.) and Misc.

Phthalate Replacements (continued)

Maleic Acid Derivatives				
Di(2-ethylhexyl)maleate (Diethyl maleate)	142-16-5	1000 µg/mL in Hexane	PLAS-PL-090S	1 mL
Di-n-butyl maleate	105-76-0	1000 µg/mL in Hexane	PLAS-PL-091S	1 mL
Mellitates				
Tricapryl trimellitate	27251-75-8	1000 µg/mL in Hexane	PLAS-PL-092S	1 mL
Triisodecyl trimellitate	36631-30-8	1000 µg/mL in Hexane	PLAS-PL-093S	1 mL
Tri(n-octyl, n-decyl) trimellitate	67989-23-5	1000 µg/mL in Hexane	PLAS-PL-094S	1 mL
Myristate				
Isopropyl myristate	110-27-0	1000 µg/mL in Hexane	PLAS-PL-095S	1 mL
Oleic Acid Derivatives				
Glycerol monooleate	25496-72-4	1000 µg/mL in Hexane	PLAS-PL-096S	1 mL
Methyl oleate	112-62-9	1000 µg/mL in Hexane	PLAS-PL-097S	1 mL
n-Propyl oleate	111-59-1	1000 µg/mL in Hexane	PLAS-PL-098S	1 mL
Tetrahydrofurfuryl oleate	5420-17-7	1000 µg/mL in Hexane	PLAS-PL-099S	1 mL
Palmitic Acid derivative				
Isopropyl palmitate	142-91-6	1000 µg/mL in Hexane	PLAS-PL-100S	1 mL
Benzoic Acid Derivatives				
Di(propylene glycol) dibenzoate	27138-31-4	1000 µg/mL in Hexane	PLAS-PL-101S	1 mL
Polyethylene glycol 200 dibenzoate	9004-86-8	1000 µg/mL in Hexane	PLAS-PL-102S	1 mL
Phosphoric Acid Derivatives				
t-Butylphenyl diphenyl phosphate	56803-37-3	1000 µg/mL in Hexane	PLAS-PL-103S	1 mL
tris(2-Butoxyethyl) phosphate	78-51-3	1000 µg/mL in Hexane	PLAS-PL-104S	1 mL
Ricinoleic Acid Derivatives				
Butyl ricinoleate	151-13-3	1000 µg/mL in Hexane	PLAS-PL-105S	1 mL
Glyceryl (triacetyl) ricinoleate	101-34-8	1000 µg/mL in Hexane	PLAS-PL-106S	1 mL
n-Butyl acetyl ricinoleate	140-04-5	1000 µg/mL in Hexane	PLAS-PL-107S	1 mL
Propylene glycol ricinoleate	26402-31-3	1000 µg/mL in Hexane	PLAS-PL-108S-T	1 mL
Succinic acid Derivatives				
Diethyl succinate	123-25-1	1000 µg/mL in Hexane	PLAS-PL-109S	1 mL
Sulfonic acid Derivatives				
o,p-Toluenesulfonamide	8013-74-9	1000 µg/mL in Hexane	PLAS-PL-110S	1 mL
N-Ethyl o,p-toluenesulfonamide	8047-99-2	1000 µg/mL in Toluene	PLAS-PL-111S-T	1 mL
Stearic acid Derivatives				
Ethylene glycol monostearate	111-60-4	1000 µg/mL in Hexane	PLAS-PL-112S	1 mL
Isopropyl isostearate	68171-33-5	1000 µg/mL in Hexane	PLAS-PL-113S	1 mL
n-Butyl stearate	123-95-5	1000 µg/mL in Hexane	PLAS-PL-114S	1 mL
Glycerol monostearate	31566-31-1	1000 µg/mL in Toluene	PLAS-PL-115S-T	1 mL
Propylene glycol monostearate	1323-39-3	1000 µg/mL in Hexane	PLAS-PL-116S	1 mL

Phthalate Set and Miscellaneous Mixtures

Phthalate Solution Set

Set of 17 Phthalate Solutions
Each at 100 µg/mL Concentration

ALR-PHT-SET 17 x 1 mL
S in MeOH, S-CN in AcCN

Compound	CAS No.	Cat. No.	Unit
Benzyl butyl phthalate	85-68-7	ALR-082S	1 mL
Diamyl phthalate	131-18-0	ALR-098S	1 mL
Dicyclohexyl phthalate	84-61-7	ALR-099S	1 mL
Di(2-ethyl hexyl) phthalate (DEHP)	117-81-7	ALR-097S	1 mL
Diethyl phthalate	84-66-2	ALR-110S	1 mL
Di-hexyl phthalate	84-75-3	ALR-100S	1 mL
Diisooctyl phthalate	26761-40-0	ALR-101S	1 mL
Diisonyl phthalate	68515-48-0	ALR-102S	1 mL
Diisooctyl phthalate	27554-26-3	ALR-103S	1 mL
Dimethyl phthalate (DMP)	131-11-3	ALR-111S	1 mL
Di-n-butyl phthalate (DBP)	84-74-2	ALR-104S	1 mL
Di-n-octyl phthalate	117-84-0	ALR-105S	1 mL
Monobenzyl phthalate (mBzP)	2528-16-7	ALR-134S-CN	1 mL
Monobutyl phthalate (mBP)	131-70-4	ALR-135S-CN	1 mL
Monoethyl phthalate (mEP)	2306-33-4	ALR-137S-CN	1 mL
Monoethylhexyl phthalate (mEHP)	4376-20-9	ALR-138S-CN	1 mL
Monomethyl phthalate	4376-18-5	ALR-139S-CN	1 mL

Appendix IX Phthalate Mix

APP-9-PHTH-MIX	1 x 1 mL
1000 µg/mL each in Cyclohexane	6 comps.
bis(2-Ethylhexyl)phthalate	
Dibutyl phthalate	
Di-n-octyl phthalate	
Diisooctyl phthalate	
Benzyl butyl phthalate	

Phthalate Mix

ASM-146	1 x 1 mL
1.0 mg/mL each in MeOH	6 comps.
Benzyl butyl phthalate	
Dimethyl phthalate	
Diethyl phthalate	
Di-n-butyl phthalate	
Di-n-octyl phthalate	
bis(2-Ethylhexyl)phthalate	

Phthalate Esters Mix

M-PHE	1 x 1 mL
M-PHE-PAK	5 x 1 mL
At stated conc.(µg/mL) in Acetone	6 comps.
Benzyl butyl phthalate	10
bis(2-Ethylhexyl)phthalate	50
Diethyl phthalate	25
Dimethyl phthalate	25
Di-n-butyl phthalate	25
Di-n-octyl phthalate	50

Phthalates EPA Methods

EPA Methods - Phthalate Standards

Method 506 Phthalate Esters by PID

Phthalate Esters

M-506	1 x 1 mL
M-506-PAK	5 x 1 mL
1.0 mg/mL each in Isooctane	7 comps.

Benzyl butyl phthalate
Dimethyl phthalate
Diethyl phthalate
Di-n-butyl phthalate

bis(2-Ethylhexyl)adipate
bis(2-Ethylhexyl)phthalate
Di-n-octyl phthalate

M-506-QC	1 x 1 mL
M-506-QC-PAK	5 x 1 mL
At stated conc. (mg/mL) in MeOH	7 comps.

Benzyl butyl phthalate	0.25	bis(2-Ethylhexyl)adipate	1.2
Dimethyl phthalate	0.1	bis(2-Ethylhexyl)phthalate	0.25
Diethyl phthalate	0.1	Di-n-octyl phthalate	0.65
Di-n-butyl phthalate	0.1		

Method 606 Phthalate Esters by GC/ECD

M-606	1 x 1 mL
M-606-PAK	5 x 1 mL

0.2 mg/mL each in MeOH

Benzyl butyl phthalate
Dimethyl phthalate
Diethyl phthalate

Di-n-butyl phthalate
Di-n-octyl phthalate
bis(2-Ethylhexyl)phthalate

Method 8060 Phthalate Esters by GC/ECD

Phthalate Esters

M-8060	1 x 1 mL
M-8060-PAK	5 x 1 mL
2.0 mg/mL each in Isooctane	6 comps.

Benzyl butyl phthalate
Diethyl phthalate
Dimethyl phthalate

Di-n-butyl phthalate
Di-n-octyl phthalate
bis(2-Ethylhexyl)phthalate

M-8060-QC	1 x 1 mL
M-8060-QC-PAK	5 x 1 mL
At stated conc. (mg/mL) in MeOH	6 comps.

Benzyl butyl phthalate	0.1	Di-n-butyl phthalate	0.25
Diethyl phthalate	0.25	Di-n-octyl phthalate	0.5
Dimethyl phthalate	0.25	bis(2-Ethylhexyl)phthalate	0.5

Method 8061A Phthalate Esters by GC/ECD

Phthalate Esters

M-8061-R1	1 x 1 mL
M-8061-R1-PAK	5 x 1 mL
1.0 mg/mL each in Hexane	15 comps.

bis(2-n-Butoxyethyl)phthalate
Butyl benzyl phthalate
Diamyl phthalate
Di-n-butyl phthalate
Dicyclohexyl phthalate
Diethyl phthalate
Dihexyl phthalate
Diisobutyl phthalate

Dimethyl phthalate
Dinonyl phthalate
Di-n-octyl phthalate
bis(2-Ethoxethyl)phthalate
bis(2-Ethylhexyl)phthalate
bis(2-Methoxethyl)phthalate
bis(4-Methyl-2-pentyl)phthalate

M-8061A	1 x 1 mL
M-8061A-PAK	5 x 1 mL
1.0 mg/mL each in Hexane	6 comps.

Butyl benzyl phthalate
bis(2-Ethylhexyl)phthalate
Di-n-butyl phthalate

Diethyl phthalate
Dimethyl phthalate
Di-n-octyl phthalate

Matrix Spike Solution

M-8061A-MS	1 x 1 mL
M-8061A-MS-PAK	5 x 1 mL
2.0 mg/mL each in Acetone	2 comps.

Butyl benzyl phthalate

bis(2-Ethylhexyl)phthalate

Internal Standard

M-8061-IS	1 x 1 mL
M-8061-IS-PAK	5 x 1 mL
5.0 mg/mL in Hexane	

Benzyl benzoate

Surrogate Standards

M-8061-SS	1 x 1 mL
M-8061-SS-PAK	5 x 1 mL

50 µg/mL each in Acetone

M-8061-SS-10X	1 x 1 mL
M-8061-SS-10X-PAK	5 x 1 mL
500 µg/mL each in Acetone	3 comps.

Dibenzyl phthalate
Diphenyl isophthalate

Diphenyl phthalate


Plastic Additives

Accelerants

Accelerators are additives that, as the name implies, accelerate or speed up the chemical reaction or the curing of the polymers into the final plastic. Accelerators are also sometimes called promoters. In rubbers, accelerators are used to increase the crosslinking reaction with sulfur in the vulcanization of rubber.

Accelerator BBTS

N-(1,1-Dimethylethyl)-2-benzothiazolesulfenamide




CAS 95-31-8 MF C₁₁H₁₄N₂S₂ MW 238.38

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-AC-003S	1 mL
NEAT	PLAS-AC-003N	50 mg

Accelerator MBT, MBT/MG

2-Mercaptobenzothiazole




CAS 149-30-4 MF C₇H₅S₂N MW 167.25

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane:Acetone 90:10	PLAS-AC-001S	1 mL
NEAT	PLAS-AC-001N	50 mg

Cure-Rite® IBT

Tetraisobutylthiuram disulfide




CAS 3064-73-1 MF C₁₈H₃₆N₂S₄ MW 408.76

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-AC-004S	1 mL
NEAT	PLAS-AC-004N	50 mg

Accelerator CBTS

N-Cyclohexyl-2-benzothioazole sulfenamide




CAS 95-33-0 MF C₁₃H₁₆N₂S₂ MW 264.41

Matrix	Cat. No.	Unit
NEAT	PLAS-AC-007N	50 mg

Accelerator EZ & EZ-SP

Zinc diethyldithiocarbamate




CAS 14324-55-1 MF C₁₀H₂₀N₂S₄ • Zn MW 361.93

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane:Acetone 90:10	PLAS-AC-006S	1 mL
NEAT	PLAS-AC-006N	50 mg

Activator OT Urea

Urea




CAS 57-13-6 MF CH₄N₂O MW 60.07

Matrix	Cat. No.	Unit
1000 µg/mL in Acetone	PLAS-AC-005S-A	1 mL
NEAT	PLAS-AC-005N	50 mg

Akroform ETU-22 PM


Ethylene thiourea



CAS 96-45-7 MF C₃H₆N₂S MW 102.11

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane:Acetone 70:30	PLAS-AC-002S	1 mL
NEAT	PLAS-AC-002N	50 mg


Dipentamethylenethiuram tetrasulfide



CAS 120-54-7 MF C₁₂H₂₀N₂S₆ MW 384.69

Matrix	Cat. No.	Unit
NEAT	PLAS-AC-009N	50 mg


1,3-Diphenyl-2-thiourea



CAS 102-08-9 MF C₁₃H₁₂N₂S MW 228.31

Matrix	Cat. No.	Unit
NEAT	PLAS-AC-008N	50 mg

1,3-Di-o-tolylguanidine



CAS 97-39-2 MF C₁₅H₁₇N₃ MW 239.32

Matrix	Cat. No.	Unit
NEAT	PLAS-AC-010N	50 mg

Property Key

CAS Chemical Abstract Service Number

MF Molecular Formula

MW Molecular Weight


Plastic Additives

Antifoams

Antifoaming agents (sometimes called defoamers) act to stop foaming during processing and typically work by reducing surface tension, breaking up the foam. Foaming can cause both processing problems and weak spots in the final product.

SF100

Dimethyl silicone fluid



CAS 9016-00-6 MF $(C_2H_6OSi)_n$ MW N/A


Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-AF-001S	1 mL
NEAT	PLAS-AF-001N	50 mg

Antidegradants

Antidegradants include a broad category of additives used in compounding to slow deterioration that can occur due to oxidation, ozone, light or any combination of these conditions. It is basically a generic term for additives that include antioxidants, antiozonants and UV Stabilizers.

Akrochem Antiox 12

Butylated reaction product of p-cresol and dicyclopentadiene




CAS 68610-51-5 MF $[C_{11}H_{20}OH \bullet C_{12}H_{23}OH]_n C_4H_9$ MW 600-800

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane:Acetone 80:20	PLAS-AD-001S	1 mL
NEAT	PLAS-AD-001N	50 mg

Ethanox® 314

1,3,5-tris(3,5-Di-tert-butyl-4-hydroxybenzyl)-1,3,5-triazine-2,4,6-(1H,3H,5H)-trione




CAS 27676-62-6 MF $C_{48}H_{69}N_3O_6$ MW 784.08

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-AX-084S	1 mL
NEAT	PLAS-AX-084N	50 mg

Santoflex® 77PD

N,N'-bis(1,4-Dimethylpentyl)-p-phenylenediamine




Flexsys

CAS 3081-14-9 MF $C_{20}H_{36}N_2$ MW 304.51

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-AD-002S	1 mL
NEAT	PLAS-AD-002N	50 mg

Santoflex® IPPD

N-Isopropyl-N'-phenyl-p-phenylenediamine




CAS 101-72-4 MF $C_{15}H_{18}N_2$ MW 226.32

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane:Acetone 80:20	PLAS-AD-003S	1 mL
NEAT	PLAS-AD-003N	50 mg

Santoflex® 6PPD

N-(1,3-Dimethylbutyl)-N'-phenyl-p-phenylenediamine



CAS 793-24-8 MF $C_{18}H_{24}N_2$ MW 268.40

Matrix	Cat. No.	Unit
NEAT	PLAS-AD-004N	50 mg

Plastic Additives

Antioxidants


Oxidation during compounding or processing can cause problems such as, loss of strength, breakdown or discoloration. Oxidation can also occur in the final product causing discolouration, scratching and loss of strength, flexibility, stiffness or gloss.

Antioxidants are used in most hydrocarbon polymers including polyethylene, polypropylene, polystyrene and ABS.

Antioxidants work to slow down the oxidation cycle, usually by scavenging free radicals. Some types of antioxidants are: organophosphites , sterically hindered phenols, amines and thioesters.

Alkanox® P27

bis(2,4-di-tert-Butylphenyl)pentaerythritol diphosphate and magnesium aluminum hydroxy carbonate hydrate




Chemtura Corporation

CAS 26741-53-7 / 11097-59-9 MF C₃₃H₅₀O₆P₂ • H₁₆Al₂Mg₆O₁₉ MW N/A

Matrix	Cat. No.	Unit
NEAT	PLAS-AX-032N	50 mg

Alkanox® TNPP

tris(Nonylphenyl) phosphite with up to 1% triisopropanol amine




CAS 26523-78-4 MF C₄₅H₆₆O₃P MW 689

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-AX-077S	1 mL
NEAT	PLAS-AX-077N	50 mg

Anox® PP18

Octadecyl 3-(3,5-di-tert-butyl-4-hydroxyphenyl) propanoate




CAS 2082-79-3 MF C₃₅H₆₂O₃ MW 530.86

Matrix	Cat. No.	Unit
NEAT	PLAS-AX-114N	50 mg

Antioxidant 60

Methyl-2-mercaptopbenzimidazole




CAS 53988-10-6 MF C₈H₈N₂S MW 164.23

Matrix	Cat. No.	Unit
1000 µg/mL in Methanol	PLAS-AX-019S-M	1 mL
NEAT	PLAS-AX-019N	50 mg

Antioxidant S

Benzenamine, N-phenyl, reaction products with 2,4,4-trimethylpentene




CAS 68411-46-1 MF C₁₂H₁₁N • C₈H₁₆ MW 393.66

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-AX-057S	1 mL
NEAT	PLAS-AX-057N	50 mg

2-(2H-Benzotriazol-2-yl)-4-(1,1,3,3-tetramethylbutyl)phenol

2-(2-Hydroxy-5-tert-octyphenyl)benzotriazole




CAS 3147-75-9 MF C₂₀H₂₅N₃O MW 323.43

Matrix	Cat. No.	Unit
NEAT	PLAS-AX-094N	50 mg

BLS® 234

2-(2-Hydroxy-3,5-di-(1,1-dimethylbenzyl))-benzotriazole




CAS 70321-86-7 MF C₃₀H₂₉N₃O MW 447.57

Matrix	Cat. No.	Unit
NEAT	PLAS-AX-088N	50 mg

BLS® 292

bis(1,2,2,6,6-pentamethyl-4-piperidinyl)sebacate and Methyl(1,2,2,6,6-pentamethyl-4-piperidinyl)sebacate




CAS 41556-26-7/8219-37-7 MF C₃₀H₅₆N₂O₄/C₂₁H₃₉NO₄ MW 508.78/369.54

Matrix	Cat. No.	Unit
NEAT	PLAS-AX-089N	50 mg

BLS® 1622

Dimethyl succinate polymer with 4-hydroxy-2,2,6,6-tetramethyl-1-piperidine ethanol



CAS 65447-77-0 MF (C₁₅H₂₅NO₄)_n MW (283.35)_n


Matrix	Cat. No.	Unit
NEAT	PLAS-AX-096N	50 mg

Plastic Additives

Antioxidants (continued)

BLS® 1944

Poly[[6-[(1,1,3,3-tetramethylbutyl)aminol]-s-triazine-2,4-diy][[(2,2,6,6-tetramethyl-4-piperidyl)imino]hexamethylene[(2,2,6,6-tetramethyl-4-piperidyl)imino]]




CAS 70624-18-9 MF (C₃₅H₆₆N₈)_n MW (599.09)_n

Matrix	Cat. No.	Unit
NEAT	PLAS-AX-090N	50 mg

BNX 1077

Benzene propanoic acid, 3,5-bis(1,1-dimethylethyl)-4-hydroxy-, isotridecyl ester




CAS 847488-62-4 MF C₃₀H₅₂O₃ MW 460.73

Matrix	Cat. No.	Unit
NEAT	PLAS-AX-087N	50 mg

BNX 1225TPR

Blend of BNX®1010, Benefos®1680 and SIS Block Copolymer




CAS 6683-19-8/31570-04-4/25038-32-8 MF N/A MW N/A

Matrix	Cat. No.	Unit
NEAT	PLAS-AX-091N	50 mg

2-tert-Butyl-6-(5-chloro-2H-benzotriazol-2-yl)-4-methylphenol

2-(2-Hydroxy-3-tert-butyl-5-methylphenyl)-5-chloro-benzotriazole




CAS 3896-11-5 MF C₁₇H₁₈ClN₃O MW 315.80

Matrix	Cat. No.	Unit
NEAT	PLAS-AX-093N	50 mg

4,4'-Butyldenebis(6-tert-butyl-m-cresol)

6,6'-di-tert-Butyl-4,4'-butyldene di-m-cresol




CAS 85-60-9 MF C₂₆H₃₈O₂ MW 382.58

Matrix	Cat. No.	Unit
NEAT	PLAS-AX-105N	50 mg

Cyanox® 1212

Lauryl stearyl thiopropionate




CAS 13103-52-1 MF C₃₈H₇₀O₄S MW 599.00

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-AX-047S	1 mL
NEAT	PLAS-AX-047N	50 mg

Cyanox® 1790

1,3,5-tris(4-tert-Butyl-3-hydroxy-2,6-dimethylbenzyl)-1,3,5-triazine-2,4,6-(1H, 3H,5H)-trione




CAS 40601-76-1 MF C₄₂H₅₇N₃O₆ MW 699.92

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-AX-005S	1 mL
NEAT	PLAS-AX-005N	50 mg

Cyanox® 2246

2,2'-Methylene-bis(4-methyl-6-tert-butylphenol)




CAS 119-47-1 MF C₂₃H₃₂O₂ MW 340.55

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-AX-013S	1 mL
NEAT	PLAS-AX-013N	50 mg

Cyanox® 425

2,2'-Methylene-bis(4-ethyl-6-tert-butylphenol)



CAS 88-24-4 MF C₂₅H₃₆O₂ MW 368.55

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-AX-012S	1 mL
NEAT	PLAS-AX-012N	50 mg

Property Key
CAS Chemical Abstract Service Number


MF Molecular Formula
MW Molecular Weight

Plastic Additives

Antioxidants (continued)

Cyanox® LTDP

Dilauryl thiodipropionate




CAS 123-28-4 MF C₃₀H₅₈O₄S MW 514.85

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-AX-041S	1 mL
NEAT	PLAS-AX-041N	50 mg

Cyanox® STDP


Distearyl thiodipropionate



CAS 693-36-7 MF C₄₂H₈₂O₄S MW 683.3

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-AX-044S	1 mL
NEAT	PLAS-AX-044N	50 mg


Dibenzylhydroxylamine



CAS 621-07-8 MF C₁₄H₁₅NO MW 213.28

Matrix	Cat. No.	Unit
NEAT	PLAS-AX-092N	50 mg

N,N'-Dibutylthiourea




CAS 109-46-6 MF C₉H₂₀N₂S MW 188.33

Matrix	Cat. No.	Unit
NEAT	PLAS-AX-102N	50 mg

3,9-bis(2,4-Dicumylphenoxy)-2,4,8,10-tetraoxa-3,9-diphosphaspiro[5,5]undecane

3,9-bis[2,4-bis(2-Phenylpropan-2-yl)phenoxy]-2,4,8,10-tetraoxa-3,9-diphosphaspiro[5,5]undecane




CAS 154862-43-8 MF C₅₃H₅₈O₆P₂ MW 852.97

Matrix	Cat. No.	Unit
NEAT	PLAS-AX-111N	50 mg

Diethyl 3,5-di-tert-butyl-4-hydroxybenzylphosphonate


2,6-di-tert-Butyl-4-(diethoxyphosphorylmethyl)phenol



CAS 976-56-7 MF C₁₉H₃₃O₄P MW 356.44

Matrix	Cat. No.	Unit
NEAT	PLAS-AX-110N	50 mg

N,N'-Diethylthiourea




CAS 105-55-5 MF C₅H₁₂N₂S MW 132.23

Matrix	Cat. No.	Unit
NEAT	PLAS-AX-103N	50 mg

3,9-bis(Octadecyloxy)-2,4,8,10-tetraoxa-3,9-diphosphaspiro[5,5]undecane

Distearyl pentaerythritol bisphosphile




CAS 3806-34-6 MF C₄₁H₈₂O₆P₂ MW 733.03

Matrix	Cat. No.	Unit
NEAT	PLAS-AX-108N	50 mg

Dipentamethylenethiuram tetrasulfide

Piperidine, 1,1'-(tetrathiodicarbonothioyl)-bis-




CAS 120-54-7 MF C₁₂H₂₀N₂S₆ MW 384.70

Matrix	Cat. No.	Unit
NEAT	PLAS-AX-104N	50 mg

1,3-Diphenyl-2-thiourea

1,3-Diphenylthiourea




CAS 102-08-9 MF C₁₃H₁₂N₂S MW 228.31

Matrix	Cat. No.	Unit
NEAT	PLAS-AX-100N	50 mg

Distyryl biphenyl

Disodium 4,4'-bis(2-sulfonatostyryl)biphenyl




CAS 27344-41-8 MF C₂₈H₂₀Na₂O₆S₂ MW 562.57

Matrix	Cat. No.	Unit
NEAT	PLAS-AX-099N	50 mg

1,3-Di-o-tolylguanidine

1,2-bis(2-Methylphenyl)guanidine




CAS 97-39-2 MF C₁₅H₁₇N₃ MW 239.32

Matrix	Cat. No.	Unit
NEAT	PLAS-AX-101N	50 mg

2,6-Di-tert-butyl-4-ethylphenol

2,6-bis(1,1-Dimethylethyl)-4-ethylphenol



CAS 4130-42-1 MF C₁₆H₂₆O MW 234.38


Matrix	Cat. No.	Unit
NEAT	PLAS-AX-107N	50 mg

Plastic Additives

Antioxidants (continued)

2,6-Di-tert-butylphenol

2,6-Di-tert-butylphenol




CAS 128-39-2 MF C₁₄H₂₂O MW 206.32

Matrix	Cat. No.	Unit
NEAT	PLAS-AX-112N	50 mg

Ethanox® 310

Pentaerythritol tetrakis (3-(3,5-di-tert-butyl-4-hydroxyphenyl)propionate)




CAS 6683-19-8 MF C₇₃H₁₀₈O₁₂ MW 1177.65

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-AX-086S	1 mL
NEAT	PLAS-AX-086N	50 mg

Ethanox® 314

tris(3,5-Di-tert-butyl-4-hydroxybenzyl)isocyanurate




CAS 27676-62-6 MF C₄₈H₆₉N₃O₆ MW 784.08

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-AX-084S	1 mL
NEAT	PLAS-AX-084N	50 mg

Ethanox® 323

Nonylphenol disulfide oligomer




CAS MF MW

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-AX-082S	1 mL
NEAT	PLAS-AX-082N	50 mg

Ethanox® 330

1,3,5-Trimethyl-2,4,6-tris(3,5-di-tert-butyl-4-hydroxybenzyl)benzene




CAS 1709-70-2 MF C₅₄H₇₈O₃ MW 775.32

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-AX-021S	1 mL
NEAT	PLAS-AX-021N	50 mg

Ethaphos® 368

Tris(2,4-di-tert-butylphenyl) phosphite




CAS 31570-04-4 MF C₄₂H₆₃O₃P MW 646.92

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-AX-074S	1 mL
NEAT	PLAS-AX-074N	50 mg

Ethanox® 376

Octadecyl 3,5-di-tert-butyl-4-hydroxyhydrocinnamate




CAS 2082-79-3 MF C₃₅H₆₂O₃ MW 530.87

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-AX-054S	1 mL
NEAT	PLAS-AX-054N	50 mg

Ethanox® 702

4,4'-Methylene bis(2,6-di-tert-butylphenol)




CAS 118-82-1 MF C₂₉H₄₄O₂ MW 424.66

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-AX-025S	1 mL
NEAT	PLAS-AX-025N	50 mg

Ethanox® 703

2,6-Di-tert-butyl-N,N-dimethylamino-p-cresol



CAS 88-27-7 MF C₁₇H₂₉NO MW 263.42

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-AX-085S	1 mL
NEAT	PLAS-AX-085N	50 mg

CAS 25973-55-1 MF C₂₂H₂₉N₃O MW 351.49


Matrix	Cat. No.	Unit
NEAT	PLAS-AX-095N	50 mg

Plastic Additives

Antioxidants (continued)

Irganox® 245

Triethyleneglycol bis[3-(3-tert-butyl-4-hydroxy-5-methylphenol)propionate]




CAS 36443-68-2 MF C₃₄H₅₀O₈ MW 586.76

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane:Acetone (95:5)	PLAS-AX-070S	1 mL
NEAT	PLAS-AX-070N	50 mg

Irganox® 259

Hexamethylene bis(3-(3,5-di-tert-butyl-4-hydroxy-phenyl)propionate)




CAS 35074-77-2 MF C₄₀H₆₂O₆ MW 638.92

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-AX-045S	1 mL
NEAT	PLAS-AX-045N	50 mg

Irganox® 565

2,4-bis(n-Octylthio)-6-(4-hydroxy-3,5-di-tert-butylanilino)-1,3,5-triazine




CAS 991-84-4 MF C₃₃H₅₆N₄OS₂ MW 588.96

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-AX-014S	1 mL
NEAT	PLAS-AX-014N	50 mg

Irganox® 1035

Thiodiethylene glycol bis(3,5-di-tert-butyl-4-hydroxyhydrocinnamate)




CAS 41484-35-9 MF C₃₈H₅₈O₆S MW 642.93

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-AX-069S	1 mL
NEAT	PLAS-AX-069N	50 mg

Irganox® 1081

6,6'-Di-tert-butyl-2,2'-thiodi-p-cresol




CAS 90-66-4 MF C₂₂H₃₀O₂S MW 358.54

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-AX-080S	1 mL
NEAT	PLAS-AX-080N	50 mg

Irganox® 1098

N,N'-1,6-Hexanediyl bis[3,5-bis(1,1-dimethylethyl)-4-hydroxy-benzenepropanamide]



CAS 23128-74-7 MF C₄₀H₆₄N₂O₄ MW 636.95

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane:Acetone 80:20	PLAS-AX-050S	1 mL
NEAT	PLAS-AX-050N	50 mg

Property Key


CAS Chemical Abstract Service Number

MF MW

Molecular Formula Molecular Weight

Irganox® 1425 WL

Ethyl 3,5-di-tert-butyl-4-hydroxybenzylphosphonate, calcium salt and polyethylene-wax mixture




CAS 65140-91-2 / 9002-88-4 MF 2C₁₇H₂₉O₄P • Ca(C₂H₄)_x MW 695

Matrix	Cat. No.	Unit
NEAT	PLAS-AX-079N	50 mg

Irganox® 3125

3,5-Di-tert-butyl-4-hydroxyhydrocinnamic triester with 1,3,5-tris[2-Hydroxyethyl]-s-triazine-2,4,6[1H,3H,5H]-trione




CAS 34137-09-2 MF C₆₀H₈₇N₃O₁₂ MW 1042.35

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane:Acetone (95:5)	PLAS-AX-020S	1 mL
NEAT	PLAS-AX-020N	50 mg

Irganox® 3114FF

1,3,5-tris(3,5-Di-tert-butyl-4-hydroxybenzyl)-1,3,5-triazine-2,4,6-(1H,3H,5H)-trione



CAS 27676-62-6 MF C₄₈H₆₉N₃O₆ MW 784.08


Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-AX-078S	1 mL
NEAT	PLAS-AX-078N	50 mg

Plastic Additives

Antioxidants (continued)

Irganox® E 201

alpha-Tocopherol




CAS 10191-41-0 MF C₂₉H₅₀O₂ MW 430.71

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-AX-027S	1 mL
NEAT	PLAS-AX-027N	50 mg

Irganox® MD 1024

1,2-bis(3,5-Di-tert-butyl-4-hydroxyhydrocinnamoyl) hydrazide




CAS 32687-78-8 MF C₃₄H₅₂N₂O₄ MW 552.79

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane:Acetone 80:20	PLAS-AX-001S	1 mL
NEAT	PLAS-AX-001N	50 mg

Isonox® 132

2,6-Di-tert-butyl-4-sec-butylphenol




CAS 17540-75-9 MF C₁₈H₃₀O MW 262.43

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-AX-018S	1 mL
NEAT	PLAS-AX-018N	50 mg

Isonox® 232

2,6-Di-tert-butyl-4-nonylphenol




CAS 4306-88-1 MF C₂₃H₄₀O MW 262.43

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-AX-063S	1 mL
NEAT	PLAS-AX-063N	50 mg

Lowinox® AH25

2,5-bis(1,1-Dimethylpropyl)-1,4-benzenediol




CAS 79-74-3 MF C₁₆H₂₆O₂ MW 250.38

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-AX-016S	1 mL
NEAT	PLAS-AX-016N	50 mg

Lowinox® CPL

Polymeric sterically hindered phenol




CAS 68610-51-5 MF N/A MW 600-700

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane:Acetone 95:5	PLAS-AX-059S	1 mL
NEAT	PLAS-AX-059N	50 mg

Lowinox® TBM-6

4,4'-Thiobis(2-tert-butyl-5-methylphenol)



CAS 96-69-5 MF C₂₂H₃₀O₂S MW 358.54

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane:Acetone 90:10	PLAS-AX-024S	1 mL
NEAT	PLAS-AX-024N	50 mg

Markstat® 60

Polyethylene glycol ether - contain < 20% sodium perchlorate


N/A

CAS 7601-89-0 MF N/A MW N/A

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-AX-028S	1 mL
NEAT	PLAS-AX-028N	50 mg

Naugard® 412S

Pentaerythritol tetrakis (beta-laurylthiopropionate)




CAS 29598-76-3 MF C₆₅H₁₂₄O₈S₄ MW 1161.94

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-AX-030S	1 mL
NEAT	PLAS-AX-030N	50 mg

Naugard® 445

4,4'-bis(alpha,alpha-Dimethylbenzyl)diphenylamine




CAS 10081-67-1 MF C₃₀H₃₁N MW 405.57

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-AX-022S	1 mL
NEAT	PLAS-AX-022N	50 mg

Naugard® 635

4-(1-Phenylethyl)-N-[4-(1-phenylethyl)phenyl]aniline



CAS 68442-68-2 MF C₃₈H₂₇N MW 377.52

Matrix	Cat. No.	Unit
NEAT	PLAS-AX-113N	50 mg

Plastic Additives

Antioxidants (continued)

Naugard® 956

Proprietary blend of primary and secondary antioxidants


N/A

CAS N/A MF N/A MW N/A

Matrix	Cat. No.	Unit
1000 µg/mL in Toluene	PLAS-AX-060S-T	1 mL
NEAT	PLAS-AX-060N	50 mg

Naugard® A

Acetone diphenylamine condensation products




CAS 68412-48-6 MF C₁₂H₁₁N • C₃H₆O MW 227.31

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane:Acetone 80:20	PLAS-AX-026S	1 mL
NEAT	PLAS-AX-026N	50 mg

Naugard® B-25

1:1 blend of Naugard 10 & Naugard 424




CAS 6683-19-8/31570-04-4 MF C₇₃H₁₀₈O₁₂• C₄₂H₆₃O₃P MW 1177.65 / 646.92

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-AX-061S	1 mL
NEAT	PLAS-AX-061N	50 mg

Naugard® BHT

2,6-Di-tert-butyl-4-methylphenol




CAS 128-37-0 MF C₁₅H₂₄O MW 220.35

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-AX-017S	1 mL
NEAT	PLAS-AX-017N	50 mg

Naugard® HM-22

Blend of phenolic primary and diphenylamine secondary antioxidants (Naugards 76 and 445)




CAS 10081-67-1/2082-79-3 MF C₃₀H₃₁N / C₃₅H₆₂O₃ MW 405.57/530.86

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-AX-033S	1 mL
NEAT	PLAS-AX-033N	50 mg

Naugard® J

N,N'-Diphenyl-p-phenylenediamine




CAS 74-31-7 MF C₁₈H₁₆N₂ MW 260.36

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane:Acetone 50:50	PLAS-AX-048S	1 mL
NEAT	PLAS-AX-048N	50 mg

Naugard® NBC

Nickel dibutyl dithiocarbamate




CAS 13927-77-0 MF C₁₈H₃₆N₂NiS₄ MW 467.45

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-AX-051S	1 mL
NEAT	PLAS-AX-051N	50 mg

Naugard® PANA

N-Phenyl-1-naphthylamine




CAS 90-30-2 MF C₁₆H₁₃N MW 219.28

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-AX-058S	1 mL
NEAT	PLAS-AX-058N	50 mg

Naugard® PHR

tris(Nonylphenyl) phosphite with up to 1% triisopropanol amine




CAS 26523-78-4 MF C₄₅H₆₉O₃P MW 689.00

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-AX-076S	1 mL
NEAT	PLAS-AX-076N	50 mg

Naugard® PS-30

Benzaniline, N-phenyl, reaction products with 2,4,4-trimethylpentene



CAS 68411-46-1 MF C₁₂H₁₁N • C₈H₁₆ MW N/A

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-AX-038S	1 mL
NEAT	PLAS-AX-038N	50 mg

Naugard® PS-35

Butylated, octylated diphenylamine-2,6 di-tert-butyl-4-sec-butyl phenol


N/A

CAS N/A MF N/A MW N/A

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-AX-046S	1 mL
NEAT	PLAS-AX-046N	50 mg

Naugard® Q Extra

1,2-Dihydro-2,2,4-trimethylquinoline (polymerized)



n=2.5-3 typically

CAS 26780-96-1 MF (C₁₂H₁₅N)_n MW (173.25)_n


Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-AX-002S	1 mL
NEAT	PLAS-AX-002N	50 mg

Plastic Additives

Antioxidants (continued)

Naugard® RM-51

tris(Nonylphenyl)phosphite




CAS 26523-78-4 MF C₄₅H₆₉C₃P MW 689.00

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-AX-034S	1 mL
NEAT	PLAS-AX-034N	50 mg

Naugard® Super Q

1,2-Dihydro-2,2,4-trimethylquinoline (polymerized)




CAS 26780-96-1 MF (C₁₂H₁₅N)_n MW (173.25)_n

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-AX-003S	1 mL
NEAT	PLAS-AX-003N	50 mg

Naugard® XL-1

2,2'-Oxamidobis[ethyl-3-(3,5-di-tert-butyl-4-hydroxy-phenyl)propionate]




CAS 70331-94-1 MF C₄₀H₆₀N₂O₈ MW 697.00

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane:Acetone 80:20	PLAS-AX-008S	1 mL
NEAT	PLAS-AX-008N	50 mg

bis(2,2,6,6-Tetramethyl-4-piperidyl) sebacate


Bis(2,2,6,6-tetramethylpiperidin-4-yl) decanedioate



CAS 52829-07-9 MF C₂₈H₅₂N₂O₄ MW 480.72

Matrix	Cat. No.	Unit
NEAT	PLAS-AX-097N	50 mg

2,2'-(2,5-Thiophenediyi)bis(5-tert-butylbenzoxazole)




CAS 7128-64-5 MF C₂₆H₂₆N₂O₂S MW 430.56

Matrix	Cat. No.	Unit
NEAT	PLAS-AX-098N	50 mg

Ultranox® 626

Bis(2,4-di-tert-butylphenyl)pentaerythritol di-phosphite




CAS 26741-53-7 MF C₃₃H₅₀O₆P₂ MW 604.62

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-AX-031S	1 mL
NEAT	PLAS-AX-031N	50 mg

Santicizer® 278

Benzyl 3-isobutyroxy-1-isopropyl-2,2-dimethylpropyl phthalate



CAS 16883-83-3 MF C₂₇H₃₄O₆ MW 454.56

Matrix	Cat. No.	Unit
NEAT	PLAS-PL-074N	50 mg

Property Key

CAS Chemical Abstract Service Number

MF Molecular Formula
MW Molecular Weight


Plastic Additives

Antiozonants

Antiozonants are materials added to plastics to slow the deterioration of the finished product that occurs from exposure to ozone. Antiozonants typically work by migrating to the surface of the product and creating an ozone-impermeable barrier or skin on the surface.

Akrochem® NIBUD

Nickel dibutyl dithiocarbamate



CAS 13927-77-0 MF C₁₈H₃₆N₂NiS₄ MW 467.45

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-AZ-001S	1 mL
NEAT	PLAS-AZ-001N	50 mg

Akrowax™ 195

A highly refined petroleum wax which is comprised of long chain saturated hydrocarbon molecules

CAS 64742-42-3 MF N/A MW N/A

Matrix	Cat. No.	Unit
NEAT	PLAS-AZ-002N	50 mg

Property Key

CAS Chemical Abstract Service Number MF Molecular Formula
MW Molecular Weight

Plastic Additives

Blowing Agents


Blowing Agents, Plasticizer

Blowing agents are sometimes also called chemical foaming agents. They are used to release gas into the plastic or resin. Blowing agents can be used to reduce weight, improve softness, provide insulation, add shock absorption properties or add resilience in the final product.

Chemical blowing agents (as opposed to physical blowing agents such as nitrogen gas) are principally organic chemicals that decompose at elevated temperatures to release a gas during decomposition that can add a cellular structure in the plastic.

CPW-100

Chlorinated paraffin wax




CAS 63449-39-8 MF Unspecified MW

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-BA-001S	1 mL
NEAT	PLAS-BA-001N	50 mg

Celogen® AZ

Carbamoyliminourea




CAS 123-77-3 MF C₂H₄N₄O₂ MW 116.08

Matrix	Cat. No.	Unit
1000 µg/mL in DMSO	PLAS-BA-002-DMSO	1 mL
NEAT	PLAS-BA-002N	50 mg

Celogen® RA

[(4-Methylphenyl)sulfonylamino]urea



CAS 10396-10-8 MF C₈H₁₁N₃O₃S MW 229.26


Matrix	Cat. No.	Unit
NEAT	PLAS-BA-003N	50 mg

Coupling Agents

Coupling agents promote the physical or chemical interaction with the polymer.

Silquest® A-187

gamma-Glycidoxypropyltrimethoxysilane




CAS 2530-83-8 MF C₉H₂₀O₅Si MW 236.38

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-CA-004S	1 mL
NEAT	PLAS-CA-004N	50 mg

Silquest® A-1102

gamma-Aminopropyltriethoxysilane (Tech grade)




CAS 919-30-2 MF C₉H₂₃NO₃Si MW 221.37

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-CA-003S	1 mL
NEAT	PLAS-CA-003N	50 mg

Silquest® A-137

Octyltriethoxysilane




CAS 2943-75-1 MF C₁₄H₃₂O₃Si MW 276.55

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-CA-005S	1 mL
NEAT	PLAS-CA-005N	50 mg

Silquest® A-1100

gamma-Aminopropyltriethoxysilane




CAS 919-30-2 MF C₉H₂₃NO₃Si MW 221.37

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-CA-002S	1 mL
NEAT	PLAS-CA-002N	50 mg

Silquest® A-1289

bis-(Triethoxysilylpropyl)tetrasulfane




CAS 40372-72-3 MF C₁₈H₄₂O₆S₄Si₂ MW 538.94

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-CA-001S	1 mL
NEAT	PLAS-CA-001N	50 mg

Silquest® A-2171

Vinylmethyldimethoxysilane



CAS 16753-62-1 MF C₅H₁₂O₂Si MW 132.24

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-CA-006S	1 mL
NEAT	PLAS-CA-006N	50 mg

Plastic Additives


Cross-Linking Agents

Crosslinking is the polymerization reaction that branches out from the main molecular chain forming a network pattern of chemical bonds. Crosslinking agents enhance this crosslinking and bonding between polymer chains.

Crosslinking adds desirable properties such as, solidity, elasticity, impermeability to gases and better electrical insulation. Crosslinking can also improve a rubber's resistance to chemicals, heat and abrasion.

F-300, F-1000, F-1500, F-2000, F-3000

Stearic acid




CAS 57-11-4 MF C₁₈H₃₆O₂ MW 284.48

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-CL-006S	1 mL
NEAT	PLAS-CL-006N	50 mg

Perkacit® MBTS

2,2'-Dithiobis(benzothiazole)




CAS 120-78-5 MF C₁₄H₈N₂S₄ MW 332.48

Matrix	Cat. No.	Unit
1000 µg/mL in Dichloromethane	PLAS-CL-001S-D	1 mL
NEAT	PLAS-CL-001N	50 mg

Perkacit® ZDEC

Zinc diethyldithiocarbamate




CAS 14324-55-1 MF C₁₀H₂₀N₂S₂Zn MW 361.9

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-CL-007S	1 mL
NEAT	PLAS-CL-007N	50 mg

Perkacit® DPG

N,N'-Diphenylguanidine




CAS 102-06-7 MF C₁₃H₁₃N₃ MW 211.27

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane:Acetone 90:10	PLAS-CL-004S	1 mL
NEAT	PLAS-CL-004N	50 mg

Perkacit® NDBC

Nickel dibutyl dithiocarbamate




CAS 13927-77-0 MF C₁₈H₃₆N₂NiS₄ MW 467.45

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-CL-005S	1 mL
NEAT	PLAS-CL-005N	50 mg

Resimene® 3520

Hexamethoxy methyl melamine




CAS 3089-11-0 MF C₁₅H₃₀N₆O₆ MW 390.51

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-CL-003S	1 mL
NEAT	PLAS-CL-003N	50 mg

Perkacit® MBT

2-Mercaptobenzothiazole



CAS 149-30-4 MF C₇H₅S₂N MW 167.25

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-CL-002S	1 mL
NEAT	PLAS-CL-002N	50 mg

Property Key

CAS Chemical Abstract Service Number

MF Molecular Formula
MW Molecular Weight

Plastic Additives


Flame Retardants

Flame retardants are added to inhibit ignition or the flammability of the end-use product. Flame retardants generally function by inhibiting the mechanisms of burning. Typical chemical elements found in compounds used as flame retardants are, aluminum, bromine, chlorine, fluorine and sulfur.

Brominated flame retardants are commonly used in polystyrene, polyesters, polyolefins, polyamides, epoxies and ABS. Decabromodiphenyl oxide is the most frequently used brominated flame retardant. The bromodiphenyl ethers are the most highly regulated of these compounds.

Some of these flame retardants are not typically added to polymers in processing, but can be found in a polymer matrix from leaching out of the contents. The largest example of this type is the Aroclors, which can often be found in a plastic matrix from having been in contact with a fluid containing these materials.


2,2',3,4,4',5',6-Heptabromodiphenyl ether



CAS 207122-16-5 MF C₁₂H₃Br₇O MW 722.48

Matrix	Cat. No.	Unit
50 µg/mL in Isooctane	BDE-183S	1 mL


2,2',4,4',5,6-Hexabromodiphenyl ether



CAS 207122-15-4 MF C₁₂H₄Br₆O MW 643.58

Matrix	Cat. No.	Unit
50 µg/mL in Isooctane	BDE-154S	1 mL


Aroclor® 1221



CAS 11104-28-2 MF Tech mix MW

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	C-221S-H-10X	1 mL
NEAT	C-221N-50MG	50 mg


2,2',4,4'-Tetrabromodiphenyl ether



CAS 5436-43-1 MF C₁₂H₆Br₄O MW 485.82

Matrix	Cat. No.	Unit
50 µg/mL in Isooctane	BDE-047S	1 mL


2,2',4,4',6-Pentabromodiphenyl ether



CAS 189084-64-8 MF C₁₂H₅Br₅O MW 564.69

Matrix	Cat. No.	Unit
50 µg/mL in Isooctane	BDE-100S	1 mL


Aroclor® 1232



CAS 11141-16-5 MF Tech mix MW

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	C-232S-H-10X	1 mL


2,2',4,4',5-Pentabromodiphenyl ether



CAS 60348-60-9 MF C₁₂H₅Br₅O MW 564.69


Matrix	Cat. No.	Unit
50 µg/mL in Isooctane	BDE-099S	1 mL

Aroclor® 1016



16% Chlorine by Mass

2,2',4,4',5,5'-Hexabromodiphenyl ether



CAS 68631-49-2 MF C₁₂H₄Br₆O MW 643.58

Matrix	Cat. No.	Unit
50 µg/mL in Isooctane	BDE-153S	1 mL


CAS 53469-21-9 MF Tech mix MW

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	C-242S-H-10X	1 mL
NEAT	C-242N-50MG	50 mg

Plastic Additives

Flame Retardants (continued)


Aroclor® 1248



CAS 12672-29-6 MF Tech mix MW N/A

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	C-248S-H-10X	1 mL
NEAT	C-248N-50MG	50 mg


Aroclor® 1262



CAS 37324-23-5 MF Tech mix MW N/A

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	C-262S-H-10X	1 mL
NEAT	C-262N-50MG	50 mg


Aroclor® 5460



CAS 11126-42-4 MF Tech mix MW N/A

Matrix	Cat. No.	Unit
35 µg/mL in Toluene	T-460S	1 mL


Aroclor® 1254



CAS 11097-69-1 MF Tech mix MW N/A

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	C-254S-H-10X	1 mL
NEAT	C-254N-50MG	50 mg


Aroclor® 1268



CAS 11100-14-4 MF Tech mix MW N/A

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	C-268S-H-10X	1 mL


Aroclor® 6050



CAS MF Tech mix MW N/A


Matrix	Cat. No.	Unit
35 µg/mL in Toluene	T-6050S	1 mL

Aroclor® 1260



CAS 11096-82-5 MF Tech mix MW N/A


Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	C-260S-H-10X	1 mL
NEAT	C-260N-50MG	50 mg



CAS 63496-31-1 MF Technical Mix MW N/A

Matrix	Cat. No.	Unit
35 µg/mL in Toluene	T-432S	1 mL


Aroclor® 5442



CAS 12642-23-8 MF Tech mix MW N/A

Matrix	Cat. No.	Unit
35 µg/mL in Toluene	T-442S	1 mL

Decabromodiphenyl ether




CAS 1163-19-5 MF C₁₂Br₁₀O MW 959.22

Matrix	Cat. No.	Unit
50 µg/mL in Isooctane	BDE-209S	1 mL

Firemaster BP4A

4,4'-(1-Methylethyldiene) bis(2,6-dibromophenol)



CAS 79-94-7 MF C₁₅H₁₂Br₄O₂ MW 543.91


Matrix	Cat. No.	Unit
100 µg/mL in Toluene	FRS-006S	1 mL

Plastic Additives

Flame Retardants (continued)

Halowax 1013

Pentachloronaphthalene




CAS 1321-64-8 MF C₁₀H₃Cl₅ MW 300.38

Matrix	Cat. No.	Unit
100 µg/mL in MeOH	N-1013S	1 mL

Halowax 1051


Octachoronaphthalene



CAS 2234-13-1 MF C₁₀Cl₈ MW 403.73

Matrix	Cat. No.	Unit
100 µg/mL in MeOH	N-1051S	1 mL

Halowax 1099




CAS 39450-05-0 MF Tech Mix MW

Matrix	Cat. No.	Unit
100 µg/mL in MeOH	N-1099S	1 mL

Saytex® 8010

Decabromodiphenyl ethane




CAS 84852-53-9 MF C₁₄H₄Br₁₀ MW 971.22

Matrix	Cat. No.	Unit
NEAT	PLAS-FR-001N	50 mg

m-Terphenyl

1,3-Diphenylbenzene




CAS 92-06-8 MF C₁₈H₁₄ MW 230.32

Matrix	Cat. No.	Unit
NEAT	T-002N	100 mg

o-Terphenyl

1,2-Diphenylbenzene




CAS 84-15-1 MF C₁₈H₁₄ MW 230.32

Matrix	Cat. No.	Unit
NEAT	T-001N	100 mg

p-Terphenyl

1,4-Diphenylbenzene



CAS 92-94-4 MF C₁₈H₁₄ MW 230.32


Matrix	Cat. No.	Unit
NEAT	T-003N	100 mg

Property Key

CAS Chemical Abstract Service Number

MF Molecular Formula
MW Molecular Weight


Tetradecachloro-m-terphenyl



CAS 42429-89-0 MF C₁₈Cl₁₄ MW 712.48

Matrix	Cat. No.	Unit
35 µg/mL in Toluene	T-005S	1 mL


Tetradecachloro-o-terphenyl



CAS MF C₁₈Cl₁₄ MW 712.48

Matrix	Cat. No.	Unit
35 µg/mL in Toluene	T-004S	1 mL

Tetradecachloro-p-terphenyl



CAS MF C₁₈Cl₁₄ MW 712.48

Matrix	Cat. No.	Unit
100 µg/mL in MeOH	T-006S	1 mL

Plastic Additives

Plasticizers


A plasticizer is a compound added to a material, usually a plastic, to make it flexible, resilient and easier to handle. Plasticizers are major components in plastics that determine the physical properties of polymer products.

Plasticizers are generally medium to high molecular weight esters of aliphatic or aromatic carboxylic acids, or sometimes of phosphoric acid. The phosphate esters are often also used for their flame retardant properties. Adipates and phthalates are also very common, but are becoming more highly regulated due to concern that they could act as endocrine disruptors.

The USEPA regulates many Phthalates and Adipates by Methods 606, 506-1 and 8061.

Benzoflex® 2-45

Diethylene glycol dibenzoate




CAS 120-55-8 MF C₁₈H₁₈O₅ MW 314.33

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-PL-015S	1 mL
NEAT	PLAS-PL-015N	50 mg

See pages 38-39 for all Bisphenol Analogs

Bisphenol A (BPA)


4,4'-Dihydroxy-2,2-diphenylpropane



CAS 80-05-7 MF C₁₅H₁₆O₂ MW 228.29

Matrix	Cat. No.	Unit
1000 µg/mL in Methanol	M-1626-01S	1 mL


n-Butyl acetyl ricinoleate



CAS 140-04-5 MF C₂₄H₄₄O₄ MW 396.60

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-PL-107S	1 mL


n-Butyl stearate



CAS 123-95-5 MF C₂₂H₄₄O₂ MW 340.58

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-PL-114S	1 mL

t-Butylphenyl diphenyl phosphate




CAS 56803-37-3 MF C₂₂H₂₃O₄P MW 382.39

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-PL-103S	1 mL

Citroflex 2


Triethyl 2-hydroxy-1,2,3-propanetricarboxylate



CAS 77-93-0 MF C₁₂H₂₀O₇ MW 276.32

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-PL-028S	1 mL
NEAT	PLAS-PL-028N	50 mg

Butyl ricinoleate




CAS 151-13-3 MF C₂₂H₄₂O₃ MW 354.57

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-PL-105S	1 mL

Celogen® SD-125

50% Azodicarbonamide in a phthalate plasticizer




CAS 77-94-1 MF C₁₈H₃₂O₇ MW 360.45

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-PL-030S	1 mL
NEAT	PLAS-PL-030N	50 mg

Citroflex A-2

Triethyl 2-acetoxy-1,2,3-propanetricarboxylate



CAS 77-89-4 MF C₁₄H₂₂O₈ MW 318.32

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-PL-001S	1 mL
NEAT	PLAS-PL-001N	50 mg


Property Key
CAS Chemical Abstract Service Number **MF** Molecular Formula
MW Molecular Weight

Plastic Additives

Plasticizers (continued)

Citroflex A-4

Tributyl 2-acetoxy-1,2,3-propanetricarboxylate




CAS 77-90-7 MF C₂₀H₃₄O₈ MW 402.54

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-PL-002S	1 mL
NEAT	PLAS-PL-002N	50 mg

Citroflex B-6


n-Butyryltri-n-hexyl citrate



CAS 82469-79-2 MF C₂₈H₅₀O₈ MW 514.7

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-PL-025S	1 mL
NEAT	PLAS-PL-025N	50 mg


Cresyl diphenyl phosphate



CAS 26444-49-5 MF C₁₉H₁₇O₄P MW 340.31

Matrix	Cat. No.	Unit
NEAT	PLAS-PL-059N	50 mg


Dibutyl phthalate



CAS 84-74-2 MF C₁₆H₂₂O₄ MW 278.34

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-PL-013S	1 mL
NEAT	PLAS-PL-013N	50 mg


Dibutyl fumarate



CAS 105-75-9 MF C₁₂H₂₀O₄ MW 228.28

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-PL-087S	1 mL


Diethyl adipate



CAS 141-28-6 MF C₁₀H₁₈O₄ MW 202.25

Matrix	Cat. No.	Unit
NEAT	PLAS-PL-043N	50 mg


Di-n-butyl maleate



CAS 105-76-0 MF C₁₂H₂₀O₄ MW 228.28

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-PL-091S	1 mL

Di(2-ethylhexyl) azelate




CAS 103-24-2 MF C₂₅H₄₈O₄ MW 412.66

Matrix	Cat. No.	Unit
1000 µg/mL in Acetone	PLAS-PL-081S-A	1 mL

Dibutyl sebacate


Dibutyl decanedioate



CAS 109-43-3 MF C₁₈H₃₄O₄ MW 314.46

Matrix	Cat. No.	Unit
NEAT	PLAS-PL-062N	50 mg


Diethyl succinate



CAS 123-25-1 MF C₈H₁₄O₄ MW 174.19

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-PL-109S	1 mL


Diisodecyl adipate



CAS 27178-16-1 MF C₂₆H₅₀O₄ MW 426.67

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-PL-083S	1 mL


Di(n-heptyl, n-nonyl) adipate



CAS 68515-75-3 MF C₂₂H₄₂O₄ MW 370.57

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-PL-080S	1 mL

Di-n-hexyl azelate




CAS 109-31-9 MF C₂₁H₄₀O₄ MW 356.54

Matrix	Cat. No.	Unit
1000 µg/mL in Acetone	PLAS-PL-078S-A	1 mL

Plastic Additives

Plasticizers (continued)


Diisobutyl adipate



CAS 141-04-8 MF C₁₄H₂₆O₄ MW 258.35

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-PL-082S	1 mL


Diisooctyl azelate



CAS 26544-17-2 MF C₂₅H₄₈O₄ MW 412.65

Matrix	Cat. No.	Unit
1000 µg/mL in Acetone	PLAS-PL-076S-A	1 mL

Diisodecyl azelate




CAS 28472-97-1 MF C₂₉H₅₆O₄ MW 468.75

Matrix	Cat. No.	Unit
1000 µg/mL in Acetone	PLAS-PL-075S-A	1 mL

Diisoctyl phthalate

bis(6-Methylheptyl)-1,2-benzenedicarboxylate




CAS 27554-26-3 MF C₂₄H₃₈O₄ MW 390.56

Matrix	Cat. No.	Unit
NEAT	PLAS-PL-071N	50 mg

Dimethyl adipate


Dimethyl hexanedioate



CAS 627-93-0 MF C₈H₁₄O₄ MW 174.19

Matrix	Cat. No.	Unit
NEAT	PLAS-PL-070N	50 mg

Dimethyl azelate




CAS 1732-10-1 MF C₁₁H₂₀O₄ MW 216.27

Matrix	Cat. No.	Unit
1000 µg/mL in Acetone	PLAS-PL-077S-A	1 mL

Dimethyl sebacate

Dimethyl decanedioate




CAS 106-79-6 MF C₁₂H₂₂O₄ MW 230.30

Matrix	Cat. No.	Unit
NEAT	PLAS-PL-061N	50 mg

Disflamoll TP

Triphenyl phosphate




CAS 115-86-6 MF C₁₈H₁₅O₄P MW 326.28

Matrix	Cat. No.	Unit
NEAT	PLAS-PL-069N	50 mg

Diethyl phthalate (DOP)

bis(2-Ethylhexyl)-1,2-benzenedicarboxylate




Houghton Chemical

CAS 117-81-7 MF C₂₄H₃₈O₄ MW 390.56

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-PL-019S	1 mL
NEAT	PLAS-PL-019N	50 mg

Di(propylene glycol) dibenzoate




CAS 27138-31-4 MF C₂₀H₂₂O₅ MW 342.39

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-PL-101S	1 mL

Disflamoll® TKP


Tricresyl phosphate



CAS 1330-78-5 MF C₂₁H₂₁O₄P MW 368.36

Matrix	Cat. No.	Unit
NEAT	PLAS-PL-073N	50 mg

Di(tridecyl) adipate



CAS 16958-92-2 MF C₃₂H₆₂O₄ MW 510.83

Matrix	Cat. No.	Unit
1000 µg/mL in Acetone	PLAS-PL-079S-A	1 mL


Epoxidized linseed oil

N/A

CAS 8016-11-3 MF N/A MW N/A

Matrix	Cat. No.	Unit
1000 µg/mL in Toluene	PLAS-PL-085S-T	1 mL

Ethylene glycol monostearate



CAS 111-60-4 MF C₂₀H₄₀O₃ MW 328.53


Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-PL-112S	1 mL

Property Key
CAS Chemical Abstract Service Number
MF Molecular Formula
MW Molecular Weight

Plastic Additives

Plasticizers (continued)

N-Ethyl-o,p-toluenesulfonamide



CAS 8047-99-2 MF C₉H₁₃NO₂S MW 199.27

Matrix	Cat. No.	Unit
1000 µg/mL in Toluene	PLAS-PL-111S-T	1 mL

2-Ethylhexyl epoxy tallate


N/A

CAS 61789-01-3 MF N/A MW N/A

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-PL-086S	1 mL

2-Ethylhexyl sebacate

bis(2-Ethylhexyl) decanedioate




CAS 122-62-3 MF C₂₆H₅₀O₄ MW 426.67

Matrix	Cat. No.	Unit
NEAT	PLAS-PL-064N	50 mg

bis(2-Ethylhexyl) terephthalate


bis(2-Ethylhexyl)-1,4-benzenedicarboxylate



CAS 6422-86-2 MF C₂₄H₃₈O₄ MW 390.56

Matrix	Cat. No.	Unit
NEAT	PLAS-PL-065N	50 mg


Glycerol monooleate



CAS 25496-72-4 MF C₂₁H₄₀O₄ MW 356.54

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-PL-096S	1 mL


Glycerol monostearate



CAS 31566-31-1 MF C₂₁H₄₂O₄ MW 358.56

Matrix	Cat. No.	Unit
1000 µg/mL in Toluene	PLAS-PL-115S-T	1 mL


Glycerol triacetate



CAS 102-76-1 MF C₉H₁₄O₆ MW 218.20

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-PL-088S	1 mL

Glyceryl (triacetyl) ricinoleate





CAS 101-34-8 MF C₆₃H₁₁₀O₁₂ MW 1059.54

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-PL-106S	1 mL

Hercoflex® 900

1,3-Isobenzofurandione, polymer with 2,2'-(1,2-ethanediylbis(oxy))bis(ethanol), benzoate



CAS 68186-30-1 MF (C₈H₄O₃)_n (C₆H₁₄O₄)_n (C₇H₆O₂)_n MW N/A

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-PL-038S	1 mL
NEAT	PLAS-PL-038N	50 mg

Hi-Point PD-1

Methyl ethyl ketone peroxide



CAS 1338-23-4 MF C₆H₁₈O MW 210.23

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-PL-024S	1 mL
NEAT	PLAS-PL-024N	50 mg


bis(2-Hydroxyethyl) dimerate

N/A

CAS 68855-78-7 MF (C₂H₆O₂₄ • unspecified)x MW N/A

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-PL-084S	1 mL


Isopropyl isostearate



CAS 68171-33-5 MF C₂₁H₄₂O₂ MW 326.56

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-PL-113S	1 mL


Isopropyl myristate



CAS 110-27-0 MF C₁₇H₃₄O₂ MW 270.45

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-PL-095S	1 mL

Isopropyl palmitate



CAS 142-91-6 MF C₁₉H₃₈O₂ MW 298.50


Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-PL-100S	1 mL

Plastic Additives

Plasticizers (continued)

Jayflex® 77

Diisoheptyl phthalate




CAS 71888-89-6 MF C₂₂H₃₄O₄ MW 362.50

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-PL-017S	1 mL
NEAT	PLAS-PL-017N	50 mg

Jayflex® L11P-E


Diundecyl phthalate



CAS 3648-20-2 MF C₃₀H₅₀O₄ MW 474.72

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-PL-021S	1 mL
NEAT	PLAS-PL-021N	50 mg


Methyl O-acetylricinoleate



CAS 140-03-4 MF C₂₁H₃₈O₄ MW 354.52

Matrix	Cat. No.	Unit
NEAT	PLAS-PL-063N	50 mg

Methyl oleate




CAS 112-62-9 MF C₁₉H₃₆O₂ MW 296.49

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-PL-097S	1 mL

Jayflex® DIDP

Diisodecyl phthalate




CAS 68515-49-1 MF C₂₈H₄₆O₄ MW 446.66

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-PL-016S	1 mL
NEAT	PLAS-PL-016N	50 mg

Jayflex® TINTM

Triisononyl trimellitate




CAS 53894-23-8 MF C₃₆H₆₀O₆ MW 588.96

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-PL-029S	1 mL
NEAT	PLAS-PL-029N	50 mg

Jayflex® DINP

Diisononyl phthalate




CAS 68515-48-0 MF C₂₆H₄₂O₄ MW 418.61

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-PL-018S	1 mL
NEAT	PLAS-PL-018N	50 mg

Jayflex® DTDP

Ditridecyl phthalate




CAS 68515-47-9 MF C₃₄H₅₈O₄ MW 530.82

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-PL-020S	1 mL
NEAT	PLAS-PL-020N	50 mg

Markstat® 51

Poly(ethylene glycol) monolaurate




CAS 9004-81-3 MF (C₂H₄O)_nC₁₂H₂₄O₂ MW N/A

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane:Toluene 90:10	PLAS-PL-003S	1 mL
NEAT	PLAS-PL-003N	50 mg

Morflex® 150

Dicyclohexyl phthalate




CAS 84-61-7 MF C₂₀H₂₆O₄ MW 330.46

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-PL-014S	1 mL
NEAT	PLAS-PL-014N	50 mg

Morflex® 190

Butylphthalyl butyl glycolate



CAS 85-70-1 MF C₁₈H₂₄O₆ MW 336.38


Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-PL-008S	1 mL
NEAT	PLAS-PL-008N	50 mg

Plastic Additives

Plasticizers (continued)

Morflex® 560

Tri-n-hexyl trimellitate




CAS 1528-49-0 MF C₂₇H₄₂O₆ MW 462.62

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-PL-031S	1 mL
NEAT	PLAS-PL-031N	50 mg

Morflex® x-1125

Triditecyl phthalate



CAS 119-06-2 MF C₃₄H₅₈O₄ MW 530.83

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-PL-033S	1 mL
NEAT	PLAS-PL-033N	50 mg

Paraplex® G-30

proprietary dibasic acid polyester mixture


N/A

CAS MF MW

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane:Toluene 88:12	PLAS-PL-027S	1 mL
NEAT	PLAS-PL-027N	50 mg

Plasthall® DINP plasticizer

Diisononyl phthalate



CAS 28553-12-0 MF C₂₆H₄₂O₄ MW 418.61

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-PL-072S	1 mL
NEAT	PLAS-PL-072N	50 mg

Plasthall® ESO


Epoxidized soybean oil

N/A

CAS 8013-07-8 MF N/A MW N/A

Matrix	Cat. No.	Unit
NEAT	PLAS-PL-035N	50 mg

n-Propyl oleate




CAS 111-59-1 MF C₂₁H₄₀O₂ MW 324.54

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-PL-098S	1 mL

Polycizer® butyl oleate

Butyl oleate




CAS 142-77-8 MF C₂₂H₄₂O₂ MW 338.57

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-PL-007S	1 mL
NEAT	PLAS-PL-007N	50 mg

Polycizer® DP 500


Dipropylene glycol dibenzoate



CAS 27138-31-4 MF C₂₀H₂₂O₅ MW 342.39

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-PL-011S	1 mL
NEAT	PLAS-PL-011N	50 mg

Polyethylene glycol 200 dibenzoate




CAS 9004-86-8 MF N/A MW N/A

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-PL-102S	1 mL

Santicizer® 141

2-Ethylhexyl diphenyl phosphate




CAS 1241-94-7 MF C₂₀H₂₇O₄P MW 362.4

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-PL-026S	1 mL
NEAT	PLAS-PL-026N	50 mg

Santicizer® 148

Mixture: isodecyldiphenyl phosphate (80-90%) / triphenyl phosphate



CAS 29761-21-5 MF C₂₂H₃₁O₄P MW 390.46


Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-PL-022S	1 mL
NEAT	PLAS-PL-022N	50 mg

Plastic Additives

Plasticizers (continued)

Santicizer® 160

Benzyl butyl phthalate




CAS 85-68-7 MF C₁₉H₂₀O₄ MW 312.37

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-PL-004S	1 mL
NEAT	PLAS-PL-004N	50 mg

Santicizer® 261

Benzyl octyl phthalate




CAS 68515-40-2 MF C₂₃H₂₈O₄ MW 368.47

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-PL-005S	1 mL
NEAT	PLAS-PL-005N	50 mg

Santicizer® 278


Benzyl 3-isobutyryloxy-1-isopropyl-2,2-dimethylpropyl phthalate



CAS 16883-83-3 MF C₂₇H₃₄O₆ MW 454.56

Matrix	Cat. No.	Unit
NEAT	PLAS-PL-074N	50 mg


Tetrahydrofurfuryl oleate



CAS 5420-17-7 MF C₂₃H₄₂O₃ MW 366.31

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-PL-099S	1 mL

o,p-Toluenesulfonamide




CAS 8013-74-9 MF C₇H₉NO₂S MW 171.22

Matrix	Cat. No.	Unit
1000 µg/mL in Acetone	PLAS-PL-110S-A	1 mL

Triisodecyl trimellitate


tris(8-Methylnonyl)-1,2,4-benzenetricarboxylate



CAS 36631-30-8 MF C₃₉H₆₆O₆ MW 630.94

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-PL-093S	1 mL

Tricapryl trimellitate




CAS 27251-75-8 MF C₃₃H₅₄O₆ MW 546.78

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-PL-092S	1 mL

Tri-butoxyethyl phosphate


tris(2-Butoxyethyl) phosphate



CAS 78-51-3 MF C₁₈H₃₉O₇P MW 398.47

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-PL-104S	1 mL


Tributyl phosphate



CAS 126-73-8 MF C₁₂H₂₇O₄P MW 266.31

Matrix	Cat. No.	Unit
NEAT	PLAS-PL-068N	50 mg

Triethyl phosphate




CAS 78-40-0 MF C₆H₁₅O₄P MW 182.15

Matrix	Cat. No.	Unit
NEAT	PLAS-PL-067N	50 mg

Plastic Additives

Plasticizers (continued)

2,2,4-Trimethyl-1,3-pentanediol-diisobutyrate




CAS 6846-50-0 MF C₁₆H₃₀O₄ MW 286.41

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-PL-089S	1 mL

2,2,4-Trimethyl-1,3-pentanediol-isobutyrate

3-Hydroxy-2,2,4-trimethylpentyl 2-methylpropanoate




CAS 25265-77-4 MF C₁₂H₂₄O₃ MW 216.32

Matrix	Cat. No.	Unit
NEAT	PLAS-PL-066N	50 mg

Trimellitate

tris(2-Ethylhexyl)-1,2,4-benzenetricarboxylic acid



CAS 3319-31-1 MF C₃₃H₅₄O₆ MW 546.78

Matrix	Cat. No.	Unit
NEAT	PLAS-PL-060N	50 mg

Vinsol® resin

gum rosin

N/A

CAS 8050-09-7

Matrix	Cat. No.	Unit
1000 µg/mL in CH ₂ Cl ₂	PLAS-PL-036S-D	1 mL
NEAT	PLAS-PL-036N	50 mg

Vinsol® powder

N/A

CAS 8050-09-7


Matrix	Cat. No.	Unit
1000 µg/mL in CH ₂ Cl ₂	PLAS-PL-037S-D	1 mL
NEAT	PLAS-PL-037N	50 mg

Property Key			
CAS	Chemical Abstract Service Number	MF	Molecular Formula
MW			Molecular Weight

Plastic Additives

Retarders (continued)


Cyanomethyl methyl(phenyl)carbamodithioate



CAS 76926-16-4 MF C₁₀H₁₀N₂S₂ MW 222.33

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-RT-009S	1 mL
NEAT	PLAS-RT-009N	50 mg


bis(Dodecylsulfanylthiocarbonyl) disulfide



CAS 870532-86-8 MF C₂₆H₅₀S₆ MW 555.07

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-RT-008S	1 mL
NEAT	PLAS-PL-008N	50 mg

bis(Thiobenzoyl) disulfide



CAS 5873-93-8 MF C₁₄H₁₀S₄ MW 306.49

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-RT-007S	1 mL
NEAT	PLAS-RT-007N	50 mg

2-(Dodecylthiocarbonothioylthio)-2-methylpropionic acid




CAS 461642-78-4 MF C₁₇H₃₂O₂S₃ MW 364.63

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-RT-010S	1 mL

Retarder AK

Phthalic anhydride



CAS 85-44-9 MF C₈H₄O₃ MW 148.12


Matrix	Cat. No.	Unit
1000 µg/mL in Hexane:Acetone(60:40)	PLAS-RT-001S	1 mL
NEAT	PLAS-RT-001N	50 mg

Stearates

Stearic acid and the metallic salts of this acid are used for many different applications depending on the polymer system. Stearates can act as lubricants, acid scavengers, anti-tack compounds, vulcanization promoters/accelerators or mold release agents.

Stearic Acid RG (rubber grade)

Stearic acid




CAS 57-11-4 MF C₁₈H₃₆O₂ MW 284.48

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-ST-001S	1 mL
NEAT	PLAS-ST-001N	50 mg

Stearic Acid TP

Stearic acid



CAS 57-11-4 MF C₁₈H₃₆O₂ MW 284.48

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-ST-002S	1 mL
NEAT	PLAS-ST-002N	50 mg

Property Key

CAS Chemical Abstract Service Number

MF Molecular Formula
MW Molecular Weight


Plastic Additives

UV Stabilizers

UV stabilizers act to protect the plastic against UV or sunlight damage such as discoloration, cracking, brittleness or other loss of desirable physical properties.

Typical UV Stabilizers are benzophenones, hindered amines and benzotriazole. Also used, but not as effective, are salicylate esters, cyanoacrylates and benzilidenes.


2-(2-Hydroxy-5-methylphenyl)benzotriazole



CAS 2440-22-4 MF C₁₃H₁₁N₃O MW 225.25

Matrix	Cat. No.	Unit
1000 µg/mL in AcCN	PLAS-UV-006S-CN	1 mL
NEAT	PLAS-UV-006N	50 mg

2-(5-tert-Butyl-2-hydroxyphenyl)benzotriazole




CAS 3147-76-0 MF C₁₆H₁₇N₃O MW 267.33

Matrix	Cat. No.	Unit
1000 µg/mL in AcCN	PLAS-UV-007S-CN	1 mL
NEAT	PLAS-UV-007N	50 mg

2-(2-Hydroxy-3,5-di(1,1-dimethylbenzyl)phenyl)benzotriazole

2-(2-Hydroxy-3,5-di-tert-amylphenyl)benzotriazole




CAS 70321-86-7 MF C₃₀H₂₉N₃O MW 447.57

Matrix	Cat. No.	Unit
1000 µg/mL in AcCN	PLAS-UV-008S-CN	1 mL
NEAT	PLAS-UV-008N	50 mg

2-tert-Butyl-6(5-chloro-2H-benzotriazol-2-yl)-4-methylphenol


2-(2-Hydroxy-3-tert-butyl-5-methylphenyl)-5-benzotriazole



CAS 3896-11-5 MF C₁₇H₁₈ClN₃O MW 315.80

Matrix	Cat. No.	Unit
1000 µg/mL in AcCN	PLAS-UV-009S-CN	1 mL
NEAT	PLAS-UV-009N	50 mg

2-(3,5-Di-tert-butyl- 2-hydroxyphenyl)benzotriazole




CAS 3846-71-7 MF C₂₀H₂₅N₃O MW 323.43

Matrix	Cat. No.	Unit
1000 µg/mL in AcCN	PLAS-UV-010S-CN	1 mL
NEAT	PLAS-UV-010N	50 mg

2,4-Di-tert-butyl-6(5-chloro-2H-benzotriazol-2-yl)phenol

2-(2-Hydroxy-3,5-di-tert-butylphenyl)-5-chlorobenzotriazole




CAS 3864-99-1 MF C₂₀H₂₄ClN₃O MW 357.88

Matrix	Cat. No.	Unit
1000 µg/mL in AcCN	PLAS-UV-011S-CN	1 mL
NEAT	PLAS-UV-011N	50 mg

2-(2H-Benzotriazol-2-yl)-4,6-di-tert-pentylphenol


2-(2-Hydroxy-3,5-di-tert-amylphenyl)benzotriazole



CAS 25973-55-1 MF C₂₂H₂₉N₃O MW 351.49

Matrix	Cat. No.	Unit
1000 µg/mL in AcCN	PLAS-UV-012S-CN	1 mL
NEAT	PLAS-UV-012N	50 mg


2-(2-Hydroxy-5-tert-octylphenyl)benzotriazole



CAS 3147-75-9 MF C₂₀H₂₅N₃O MW 323.43

Matrix	Cat. No.	Unit
1000 µg/mL in AcCN	PLAS-UV-013S-CN	1 mL
NEAT	PLAS-UV-013N	50 mg

2-(3-sec-Butyl-5-tert-butyl-2-hydroxy-phenyl)benzotriazole




CAS 36437-37-3 MF C₂₀H₂₅N₃O MW 323.43

Matrix	Cat. No.	Unit
1000 µg/mL in AcCN	PLAS-UV-014S-CN	1 mL
NEAT	PLAS-UV-014N	50 mg

2-(2H-Benzotriazol-2-yl)-4-methyl,-6-(2-propenyl)phenol

2-(2-Hydroxy-3-allyl-5-methylphenyl)benzotriazole



CAS 2170-39-0 MF C₁₆H₁₅N₃O MW 265.31

Matrix	Cat. No.	Unit
1000 µg/mL in AcCN	PLAS-UV-015S-CN	1 mL
NEAT	PLAS-UV-015N	50 mg


UV Stabilizer Set (Solutions)	10 x 1 mL
PLAS-UV-STAB-SET	PLAS-UV -006S-CN to 015S-CN

Plastic Additives

UV Stabilizers (continued)

Tinuvin® PED

2-(2-Hydroxy-5-methylphenyl)benzotriazole




CAS 2440-22-4 MF C₁₃H₁₁N₃O MW 225.27

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-UV-005S	1 mL
NEAT	PLAS-UV-005N	50 mg

Uvinul® 3008

2-Hydroxy-4-octyloxybenzophenone




CAS 1843-05-6 MF C₂₁H₂₆O₃ MW 326.43

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-UV-002S	1 mL
NEAT	PLAS-UV-002N	50 mg

Uvinul® 3049

2,2'-Dihydroxy-4,4'-dimethoxybenzophenone




CAS 131-54-4 MF C₁₅H₁₄O₅ MW 274

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane:Acetone (80:20)	PLAS-UV-004S	1 mL
NEAT	PLAS-UV-004N	50 mg

Uvinul® 3000

2,4-Dihydroxybenzophenone




CAS 131-56-6 MF C₁₃H₁₀O₃ MW 214.22

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane:Acetone (80:20)	PLAS-UV-001S	1 mL
NEAT	PLAS-UV-001N	50 mg

Uvinul® 3040

2-Hydroxy-4-methoxybenzophenone



CAS 131-57-7 MF C₁₄H₁₂O₃ MW 228.26

Matrix	Cat. No.	Unit
1000 µg/mL in Hexane	PLAS-UV-003S	1 mL
NEAT	PLAS-UV-003N	50 mg

Vegetable Oils

Vegetable oils, typically the epoxide or the ester of the parent oil, are used as plasticizers. They not only offer the advantage of providing flexibility in the final plastic, but also add heat and light stabilizing advantages without the requirements for additional additives. Vegetable oil plasticizers are generally less toxic than their petrochemical counterparts, which makes them very attractive for certain applications like food or toys.

Some of their disadvantages are that they may not mix properly at higher concentrations, may cause brittleness in some applications and often are only suitable as secondary plasticizers.

Akrofax™ A

Vulcanized vegetable oil

N/A

CAS 68952-47-6 MF N/A MW N/A

Matrix	Cat. No.	Unit
NEAT	PLAS-VA-001N	50 mg

Akrofax™ B

Vulcanized vegetable oil

N/A

CAS N/A MF N/A MW N/A

Matrix	Cat. No.	Unit
NEAT	PLAS-VA-002N	50 mg

Property Key

CAS	Chemical Abstract Service Number	MF	Molecular Formula
		MW	Molecular Weight

Bisphenol Analogs

Bisphenol Analog

Bisphenol A (2,2'-bis(4-Hydroxyphenyl)propane, BPA) has been used in commercial and industrial applications since the 1970's. It has been the subject of numerous toxicological studies due to human exposure from leachate originating from polycarbonate plastics and epoxy-lined food and drink containers.

The evidence of the toxic effects of BPA has led to restrictions and regulations, resulting in its replacement in commercial products with related compounds. Several chemicals with structural similarity to BPA (i.e. two hydroxyl phenyl moieties) have been used as alternatives in the manufacture of polycarbonate plastics and epoxy resins. 4,4-Sulfonyldiphenol (BPS) and 4,4-Dihydroxydiphenylmethane (BPF) are the two main substitutes. However, their similarity to BPA has led to their monitoring and testing for human exposure and toxicity as well.

In addition to the BPA analogs, there has been increased scrutiny of bisphenol A diglycidyl ether (BADGE) which is a widely used building block of epoxy resin. Studies have shown that it also might be linked to adverse human health effects.


These BPA related compounds of reference standards for BPA analogs as well as the BADGE starting material is a comprehensive listing.

References:

1. Environ. Sci. Technol. 2012, 46, 9138-9145
2. Environ. Sci. Technol. 2012, 46, 12968-12976
3. Environ. Sci. Technol. 2012, 46, 11558-11565

Bisphenol A diglycidyl ether (BADGE)

2,2-bis(4-Glycidyloxyphenyl)propane




CAS 1675-54-3 MF C₂₁H₂₄O₄ MW 340.41
log Kow 3.8 PS L SOL MeOH SG 1.16 g/cm³
MP N/A BP >200 °C FP N/A

Matrix	Cat. No.	Unit
NEAT	BADGE-001N	50 mg
10 mg/mL in MeOH	BADGE-001S	1 mL

Bisphenol A (BPA)

2,2-bis(4-Hydroxyphenyl)propane




CAS 80-05-7 MF C₁₅H₁₆O₂ MW 228.29
log Kow 3.43 PS S SOL Acetone, Benzene, Ether
SG 1.14 g/cm³ MP 156 °C BP 220 °C FP 192 °C

Matrix	Cat. No.	Unit
NEAT	BPA-A-N	50 mg
10 mg/mL in MeOH	BPA-A-S	1 mL

Bisphenol AF

2,2-bis(4-Hydroxyphenyl)-1,1,1,3,3,3-hexafluoropropane




CAS 1478-61-1 MF C₁₅H₁₀F₆O₂ MW 336.23
log Kow 4.47 PS S SG 1.45 g/cm³ MP 159-162 °C
BP 400 °C FP 162 °C

Matrix	Cat. No.	Unit
NEAT	BPA-AF-N	50 mg
10 mg/mL in MeOH	BPA-AF-S	1 mL

Bisphenol AP

1,1-bis(4-Hydroxyphenyl)-1-phenylethane




CAS 1571-75-1 MF C₂₀H₁₈O₂ MW 290.36
log Kow 4.86 PS S SG 1.18 g/cm³ MP 182-183 °C
BP 473-475 °C FP 222 °C

Matrix	Cat. No.	Unit
NEAT	BPA-AP-N	50 mg
10 mg/mL in MeOH	BPA-AP-S	1 mL

Bisphenol BP

bis(4-Hydroxyphenyl)diphenylmethane




CAS 1844-01-5 MF C₂₅H₂₀O₂ MW 352.43
log Kow 6.08 PS S SG 1.20 g/cm³ MP 216-217 °C
BP 508-510 °C FP 241 °C

Matrix	Cat. No.	Unit
NEAT	BPA-BP-N	50 mg
10 mg/mL in MeOH	BPA-BP-S	1 mL

Bisphenol B

2,2-bis(4-Hydroxyphenyl)butane




CAS 77-40-7 MF C₁₆H₁₈O₂ MW 242.31
log Kow 4.13 PS S SG 1.12 g/cm³ MP 126 °C
BP 412-414 °C FP 196 °C

Matrix	Cat. No.	Unit
NEAT	BPA-B-N-10MG	10 mg
10 mg/mL in MeOH	BPA-B-S	1 mL

Bisphenol C

2,2-bis(4-Hydroxy-3-methylphenyl)propane



CAS 79-97-0 MF C₁₇H₂₀O₂ MW 256.34
log Kow 4.74 PS S SG 1.15 g/cm³ MP 139 °C
BP 390 °C FP >190 °C

Matrix	Cat. No.	Unit
NEAT	BPA-C-N	50 mg
10 mg/mL in MeOH	BPA-C-S	1 mL


Property Key

CAS	Chemical Abstract Service Number
MF	Molecular Formula
MW	Molecular Weight
PS	Physical State (Solid, Liquid)
SOL	Solubility
SG	Specific Gravity (g/cm ³)
MP	Melting Point (°C)
BP	Boiling Point (°C)
FP	Flash Point (°C)

Bisphenol Analogs

Bisphenol Analog (continued)


Bisphenol C-dichloride



CAS 14868-03-2 MF $C_{14}H_{10}Cl_2O_2$ MW 281.13
log Kow 3.75 PS S SG 1.45 g/cm³ MP 216 °C
BP 395-398 °C FP >200 °C

Matrix	Cat. No.	Unit
NEAT	BPA-C2-N	20 mg
10 mg/mL in MeOH	BPA-C2-S	1 mL


Bisphenol E



CAS 2081-08-5 MF $C_{14}H_{14}O_2$ MW 214.26
log Kow 3.19 PS S SG 1.25 g/cm³ MP 125 °C
BP 350-370 °C FP >180 °C

Matrix	Cat. No.	Unit
NEAT	BPA-E-N	50 mg
10 mg/mL in MeOH	BPA-E-S	1 mL


Bisphenol F



CAS 620-92-8 MF $C_{13}H_{12}O_2$ MW 200.23
log Kow 2.91 PS S SG 1.21 g/cm³ MP 163 °C
BP 389-390 °C FP 193 °C

Matrix	Cat. No.	Unit
NEAT	BPA-F-N-10MG	10 mg
10 mg/mL in MeOH	BPA-F-S	1 mL


Bisphenol G



CAS 127-54-8 MF $C_{21}H_{28}O_2$ MW 312.45
log Kow 6.55 PS S SG 1.05 g/cm³ MP 98 °C
BP 419-420 °C FP 185 °C

Matrix	Cat. No.	Unit
NEAT	BPA-G-N	20 mg
10 mg/mL in MeOH	BPA-G-S	1 mL


Bisphenol M



CAS 13595-25-0 MF $C_{24}H_{26}O_2$ MW 346.46
log Kow 6.25 PS S SG 1.15 g/cm³ MP 138 °C
BP >495 °C FP >200 °C

Matrix	Cat. No.	Unit
NEAT	BPA-M-N	20 mg
10 mg/mL in MeOH	BPA-M-S	1 mL

Bisphenol S




CAS 80-09-1 MF $C_{12}H_{10}O_4S$ MW 250.27
log Kow 1.65 PS S SG 1.43 g/cm³ MP 245-250 °C
BP 505-506 °C FP 259 °C

Matrix	Cat. No.	Unit
NEAT	BPA-S-N	50 mg
10 mg/mL in MeOH	BPA-S-S	1 mL

Bisphenol TMC


1,1-bis(4-Hydroxyphenyl)-3,3,5-trimethylcyclohexane



CAS 129188-99-4 MF $C_{21}H_{26}O_2$ MW 310.43
log Kow 6.29 PS S SG 1.10 g/cm³ MP 204-207 °C
BP 450 °C FP 203 °C

Matrix	Cat. No.	Unit
NEAT	BPA-TMC-N-10MG	10 mg
10 mg/mL in MeOH	BPA-TMC-S	1 mL


Bisphenol PH



CAS 24038-68-4 MF $C_{24}H_{24}O_2$ MW 380.48
log Kow 7.17 PS S SG 1.20 g/cm³ MP 118 °C
BP 567-568 °C FP 250 °C

Matrix	Cat. No.	Unit
NEAT	BPA-PH-N	20 mg
10 mg/mL in MeOH	BPA-PH-S	1 mL

Bisphenol Z



CAS 843-55-0 MF $C_{18}H_{20}O_2$ MW 268.35
log Kow 5.00 PS S SG 1.17 g/cm³ MP 189-192 °C
BP 440-441 °C FP 207 °C

Matrix	Cat. No.	Unit
NEAT	BPA-Z-N	50 mg
10 mg/mL in MeOH	BPA-Z-S	1 mL

Plastic Packaging Testing

ASTM Method D6042-92 Plastic Packaging Testing Standards

This method is used by both pharmaceutical companies and plastics manufacturers to ensure the quality of the plastic product during the manufacturing process. Compounds are often added to the method's analyte list by pharmaceutical companies.

Calibration Mix

PLAS-CAL-001	1 x 1 mL
PLAS-CAL-001-PAK	SAVE 5 x 1 mL
50 µg/mL each in Isopropanol	7 comps.
BHT	Irganox 3114
Erucamide Slip	Irganox 1010
Vitamin E	Irganox 1076
Irgafos 168	

Internal Standard Mix


PLAS-IS-001	1 x 1 mL
PLAS-IS-001-PAK	SAVE 5 x 1 mL
51.8 µg/mL in Isopropanol	

Tinuvin P

Expanded List of Additives

Each at 50 µg/mL in Isopropanol


Ultranox 626	PLAS-CAL-002-1	1 mL
Santanol R	PLAS-CAL-002-2	1 mL
Ethanox 330	PLAS-CAL-002-3	1 mL
Ethanox 323	PLAS-CAL-002-4	1 mL
Ethanox 702	PLAS-CAL-002-5	1 mL
Ethanox 703	PLAS-CAL-002-6	1 mL
Irganox 1035	PLAS-CAL-002-7	1 mL



The figure shows the separation of the compounds on the method's analyte list, as analyzed by our HPLC specialists. The primary calibration standard mixture contains the common antioxidants and slips listed in ASTM D6042-96.

The perfect companion for your analysis!

This reference book contains the compounds in this catalog, with important reference data to aid in testing and compliance.



Each Compound has:

Chemical Information

- Structure
- CAS Number (where applicable)
- RTECS Number (where available)
- Formula
- Molecular Weight
- IUPAC Name, other common names and some popular brand names

Other Important Information

- Application
- Regulatory
- Environmental Impact
- Point of Release
- Toxicological Data

Physical Properties

- Appearance
- Melting and Boiling Points
- Stability
- Solubilities in several common solvents

Analytical Data

- Mass Spectrum with key ions tabulated
- Chromatogram with conditions

Cat. No: BOOK-PLAS-002

Handbook for the Chemical Analysis of Plastic and Polymer Additives, 2nd Ed.

The Second Edition of this handbook provides the necessary tools for chemists to obtain a more complete listing of additives present in a particular polymeric matrix. This edition features:

- Updated material to include the most recent additives available
- Actual analytical data for each chemical along with the description and methods used for obtaining the results
- Toxicological and environmental impact of each product highlighted
- Summaries of regulatory and health information in a convenient "one-step" format

With 50 additional compounds, this 2nd edition nearly doubles the number of additives in several categories including processing aids, anti-static compounds, mold release products and blowing agents. It includes a listing that can be cross-referenced by trade name, chemical name, CAS number and even key mass unit ions from the GC/MS run. Also included are case studies related to "real-world" issues, tips for analysis in challenging matrices and more.

Dyes and Breakdown Products

Dyes and colorant products are one of the largest categories of plastic additives and are also used in textiles, leather goods, food and personal care products. They are used for both aesthetic purposes and to alter physical properties of the product, e.i. repelling light. Many dyes and their breakdown products have been determined to have both adverse health properties and adverse environmental properties, and as such, are being increasingly regulated. EU Directives 67/548/EEC and 2002/61/EC and 76/768/EEC are the most far-reaching regulations for this class of compounds.

EU Directive 67/548/EEC Dyes

Criterion #22 Regulated Dyes - Carcinogenic Each in 100 µg/mL in MeOH		
	Cat. No.	Unit
Disperse Blue 1	DYE-001S	1 mL
Disperse Orange 11	DYE-002S	1 mL
Disperse Yellow 3	DYE-003S	1 mL
Basic Violet 14	DYE-012S	1 mL
Direct Black 38	DYE-013S	1 mL
Direct Blue 6	DYE-014S	1 mL

Criterion #23 Regulated Dye - Disperse dyes, Sensitizing Each in 100 µg/mL in MeOH		
	Cat. No.	Unit
Disperse Blue 3	DYE-004S	1 mL
Disperse Orange 1	DYE-005S	1 mL
Disperse Orange 3	DYE-006S	1 mL
Disperse Red 1	DYE-007S	1 mL
Disperse Yellow 9	DYE-008S	1 mL
Disperse Blue 35	DYE-009S	1 mL
Disperse Blue 124	DYE-010S	1 mL
Disperse Orange 37	DYE-011S	1 mL
Disperse Blue 7	DYE-015S	1 mL
Disperse Blue 26	DYE-016S	1 mL
Disperse Blue 102	DYE-017S	1 mL
Disperse Red 11	DYE-018S	1 mL
Disperse Red 17	DYE-019S	1 mL

EU Directive 2002/61/EC Determination of Aryl Amine Breakdown Products in Azo Dyes

Individual Aryl Amine Standards

Compound	100 µg/mL in AcCN 1 mL	1000 µg/mL in AcCN 1 mL	10 µg/mL in Ethyl acetate for 10 mL 10 mL
	RAC-01	RAC-01-10X	RAC-01-EA-0.1X-10ML
o-Aminoazotoluene	RAC-02	RAC-02-10X	RAC-02-EA-0.1X-10ML
4-Aminobiphenyl	RAC-03	RAC-03-10X	RAC-03-EA-0.1X-10ML
2-Amino-4-nitrotoluene	RAC-04	RAC-04-10X	RAC-04-EA-0.1X-10ML
Benzidine †	RAC-05	RAC-05-10X	RAC-05-EA-0.1X-10ML
4-Chloroaniline	RAC-06	RAC-06-10X	RAC-06-EA-0.1X-10ML
4-Chloro-o-toluidine	RAC-07	RAC-07-10X	RAC-07-EA-0.1X-10ML
p-Cresidine	RAC-08	RAC-08-10X	RAC-08-EA-0.1X-10ML
2,4-Diaminoanisole *	RAC-09	RAC-09-10X	RAC-09-EA-0.1X-10ML
4,4'-Diaminodiphenylmethane	RAC-10	RAC-10-10X	RAC-10-EA-0.1X-10ML
2,4-Diaminotoluene	RAC-11	RAC-11-10X	RAC-11-EA-0.1X-10ML
3,3'-Dichlorobenzidine †	RAC-12	RAC-12-10X	RAC-12-EA-0.1X-10ML
3,3'-Dimethoxybenzidine †	RAC-13	RAC-13-10X	RAC-13-EA-0.1X-10ML
3,3'-Dimethylbenzidine †	RAC-14	RAC-14-10X	RAC-14-EA-0.1X-10ML
3,3'-Dimethyl-4,4'-diaminodiphenylmethane	RAC-15	RAC-15-10X	RAC-15-EA-0.1X-10ML
4,4'-Methylenebis(2-chloroaniline)	RAC-16	RAC-16-10X	RAC-16-EA-0.1X-10ML
2-Naphthylamine	RAC-17	RAC-17-10X	RAC-17-EA-0.1X-10ML
4,4'-Oxydianiline	RAC-18	RAC-18-10X	RAC-18-EA-0.1X-10ML
4,4'-Thiodianiline	RAC-19	RAC-19-10X	RAC-19-EA-0.1X-10ML
o-Toluidine	RAC-20	RAC-20-10X	RAC-20-EA-0.1X-10ML
2,4,5-Trimethylaniline	RAC-21	RAC-21-10X	RAC-21-EA-0.1X-10ML
p-Aminoazobenzene	RAC-22	RAC-22-10X	RAC-22-EA-0.1X-10ML
2-Aminobiphenyl	RAC-23	RAC-23-10X	RAC-23-EA-0.1X-10ML
o-Anisidine	RAC-24	RAC-24-10X	RAC-24-EA-0.1X-10ML

RAC-R1-SET

24 x 1 mL

100 µg/mL

* In the form of the Sulfate hydrate 171 µg/mL in Pyridine (100 µg/mL as the base)

RAC-R1-10X-SET

24 x 1 mL

1000 µg/mL

* In the form of the Sulfate hydrate 1,710 µg/mL in Pyridine (1000 µg/mL as the base)

† Subject to oxidation

Carcinogenic Aryl Amine Mix

AE-00049-R1
10 µg/mL in Ethyl acetate
23 comps.

AE-00049-R1-10ML
10 µg/mL in Ethyl acetate
23 comps.

o-Aminoazotoluene
4-Aminobiphenyl
2-Amino-4-nitrotoluene
Benzidine †
4-Chloroaniline
4-Chloro-o-toluidine
p-Cresidine
4,4'-Diaminodiphenylmethane
2,4-Diaminotoluene
3,3'-Dichlorobenzidine †
3,3'-Dimethoxybenzidine †
3,3'-Dimethylbenzidine †
3,3'-Dimethyl-4,4'-diaminodiphenylmethane
4,4'-Methylenebis(2-chloroaniline)
2-Naphthylamine
4,4'-Oxydianiline
4,4'-Thiodianiline
o-Toluidine
2,4,5-Trimethylaniline
p-Aminoazobenzene
2-Aminobiphenyl
o-Anisidine
3-Chloro-o-toluidine

Internal Standards

RAC-IS
1000 µg/mL in AcCN

RAC-IS-EA
1000 µg/mL in Ethyl acetate

3,3',5,5'-Tetramethylbenzidine †

AE-00049-SET
2 x 1 mL
AE-00049-R1, RAC-08

Plastic Additive Index

A

Accelerator BBTS 10
Accelerator CBTs 10
Accelerator EZ & EZ-SP 10
Accelerator MBT, MBT/MG 10
Activator OT Urea 10
Akrochem Antiox 12 11
Akrochem® Ceresin Wax 34
Akrochem® NIBUD 20
Akrochem® Retarder BAX 34
Akrofax™ A 37
Akrofax™ B 37
Akroform ETU-22 PM 10
Akrowax™ 195 20
Alkanox® P27 12
Alkanox® TNPP 12
2-Amino-4-nitrotoluene 41
o-Aminoazotoluene 41
p-Aminoazobenzene 41
2-Aminobiphenyl 41
4-Aminobiphenyl 41
o-Anisidine 41
Anox® PP18 12
Antioxidant 60 12
Antioxidant S 12
Aroclor® 1016 23
Aroclor® 1221 23
Aroclor® 1232 23
Aroclor® 1242 23
Aroclor® 1248 24
Aroclor® 1254 24
Aroclor® 1260 24
Aroclor® 1262 24
Aroclor® 1268 24
Aroclor® 5432 24
Aroclor® 5442 24
Aroclor® 5460 24
Aroclor® 6050 24

B

Basic Violet 14 41
Benzidine 41
Benzoflex® 2-45 26
2-(2H-Benzotriazol-2-yl)-4-(1,1,3,3-tetramethylbutyl)phenol 12
Benzyl 2-ethylhexyl phthalate 7
Benzyl benzoate 9
Benzyl butyl phthalate 2, 8, 9
bis(2,2,6,6-Tetramethyl-4-piperidyl) sebacate 19
bis(2-Hydroxyethyl) dimerate 7, 29
bis(2-Methoxyethyl)phthalate 2, 9
bis(2-n-Butoxyethyl)phthalate 2, 9
bis(4-Methyl-2-pentyl)phthalate 2, 9
bis(Dodecylsulfanylthiocarbonyl) disulfide 35
Bisphenol A (BPA) 26, 38
Bisphenol A diglycidyl ether (BADGE) 38
Bisphenol AF 38
Bisphenol AP 38
Bisphenol B 38
Bisphenol BP 38
Bisphenol C 38
Bisphenol C-dichloride 39
Bisphenol E 39
Bisphenol F 39
Bisphenol G 39
Bisphenol M 39
Bisphenol P 39
Bisphenol PH 39
Bisphenol S 39
Bisphenol TMC 39
Bisphenol Z 39
BLS® 234 12
BLS® 292 12
BLS® 1622 12
BLS® 1944 13
BNX 1077 13
BNX 1225 13

n-Butyl acetyl ricinoleate 8, 26

n-Butyl benzyl phthalate 7
Butyl benzyl phthalate 9
Butyl cyclohexyl phthalate 7
n-Butyl isobutyl phthalate 7
Butyl octyl phthalate 7
t-Butylphenyl diphenyl phosphate 8
Butyl ricinoleate 8, 26
n-Butyl stearate 8, 26
2-tert-Butyl-6-(5-chloro-2H-benzotriazol-2-yl)-4-methylphenol 13
4,4'-Butyldienebis(6-tert-butyl-m-cresol) 13

C

Celogen® AZ 21
Celogen® RA 21
Celogen® SD-125 26
4-Chloroaniline 41
3-Chloro-o-toluidine 41
4-Chloro-o-toluidine 41
Citroflex 2 26
Citroflex 4 26
Citroflex A-2 26
Citroflex A-4 27
Citroflex B-6 27
CPW-100 21
p-Cresidine 41
Cresyl diphenyl phosphate 27
Cure-Rite® IBT 10
4-Cyano-4-[(dodecylsulfanylthiocarbonyl)sulfanyl]pentanoic acid 34
4-Cyano-4-(phenylcarbonothioylthio)pentanoic acid 34
2-Cyano-2-propyl benzodithioate 34
2-Cyano-2-propyl dodecyl trithiocarbonate 34
Cyanomethyl dodecyl trithiocarbonate 34
Cyanomethyl methyl(phenyl)carbamodithioate 35
Cyanox® 425 13
Cyanox® 1212 13
Cyanox® 1790 13
Cyanox® 2246 13
Cyanox® LTDP 14
Cyanox® STDP 14

D

Decabromodiphenyl ether 24
Decyl octyl phthalate 7
Di(2-ethyl hexyl) azelate 7
Di(2-ethylhexyl) maleate 27
Di(2-ethylhexyl)maleate 8
Di(2-ethyl hexyl) phthalate (DEHP) 8
2,6-Di-tert-butyl-4-ethylphenol 14
2,6-Di-tert-butylphenol 15
Diallyl phthalate 2
4,4'-Diaminodiphenylmethane 41
2,4-Diaminoanisole 41
2,4-Diaminotoluene 41
Diamyl phthalate 2, 8, 9
Dibenzylhydroxylamine 14
Dibenzyl phthalate 2, 9
Dibutyl fumarate 7, 27
Dibutyl phthalate 2, 8, 27
Dibutyl sebacate 27
3,3'-Dichlorobenzidine 41
Dicyclohexyl phthalate 2, 8, 9
Dicyclohexyl phthalate-3,4,5,6-d4 6
Didecyl phthalate (Tech Mix) 7
Didodecyl phthalate 3
Diethyl 3,5-di-tert-butyl-4-hydroxybenzylphosphonate 14
Diethyl adipate 27
Diethyl phthalate 3, 8, 9
Diethyl phthalate-3,4,5,6-d4 6
Diethyl succinate 8, 27
Diethyl terephthalate 4
Di-hexyl phthalate 8
Dihexyl phthalate 3, 9

Diisobutyl adipate 7, 28

Diisobutyl phthalate 3, 9
Diisobutyl phthalate-3,4,5,6-d4 6
Diisooctyl azelate 7, 28
Diisodecyl adipate 7, 27
Diisodecyl azelate 7, 28
Diisodecyl phthalate 7, 8
Diisooheptyl phthalate 7
Diisohexyl phthalate 7
Diisonyl phthalate 8
Diisonyl phthalate (C8 to C10 Isomers) 7
Diisoctyl phthalate 8, 28
Diisoctyl phthalate (C8 Isomers) 7
Diisonyl phthalate 8

Diisopentyl phthalate 3
Diisopropyl phthalate 3
Dimethyl adipate 28
Dimethyl azelate 7, 28
Dimethyl isophthalate 4
Dimethyl phthalate 3, 8, 9
Dimethyl phthalate-3,4,5,6-d4 6
Dimethyl phthalate (DMP) 8
Dimethyl sebacate 28
Dimethyl terephthalate 4
Di-n-butyl maleate 8, 27
Di-n-butyl phthalate 8, 9
Di-n-butyl phthalate-3,4,5,6-d4 6
Di-n-butyl phthalate (DBP) 8
Di(n-heptyl, n-nonyl) adipate 7, 27
Di-n-heptyl phthalate 3
Di-n-hexyl azelate 7, 27
Di-n-hexyl phthalate-3,4,5,6-d4 6
Di-n-octyl phthalate 3, 8, 9
Di-n-octyl phthalate-3,4,5,6-d4 6
Dinonyl phthalate 7, 9
Di-n-pentyl phthalate-3,4,5,6-d4 6
Di-n-propyl phthalate-3,4,5,6-d4 6
Diocetyl phthalate 28
1,3-Di-o-tolylguanidine 10, 14
3,3'-Dimethoxybenzidine 41
3,3'-Dimethyl-4,4'-diaminodiphenylmethane 41
3,3'-Dimethylbenzidine 41
Dipentamethylenethiuram 14
Dipentamethylenethiuram tetrasulfide 10
Diphenyl isophthalate 4, 9
Diphenyl phthalate 3, 9
1,3-Diphenyl-2-thiourea 10, 14
Di(propylene glycol) dibenzoate 8, 28

Direct Black 38 41

Disflamoll® TKP 28

Disflamoll TP 28

Disperse Blue 1 41

Disperse Blue 3 41

Disperse Blue 7 41

Disperse Blue 26 41

Disperse Blue 35 41

Disperse Blue 102 41

Disperse Blue 124 41

Disperse Orange 1 41

Disperse Orange 3 41

Disperse Orange 11 41

Disperse Orange 37 41

Disperse Red 1 41

Disperse Red 11 41

Disperse Red 17 41

Disperse Yellow 3 41

Disperse Yellow 9 41

Distyryl biphenyl 14

Di(tridecyl) adipate 7, 28

Diundecyl phthalate 3

2-(Dodecylthiocarbonothioylthio)-2-methylpropionic acid 35

Epoxidized linseed oil 7, 28

Ethanox® 310 15

Ethanox® 314 11, 15

E

Plastic Additive Index

Ethanox® 323 15, 40

Ethanox® 330 15, 40

Ethanox® 376 15

Ethanox® 702 15, 40

Ethanox® 703 15, 40

Ethaphos® 368 15

bis(2-Ethoxyethyl) phthalate 9

Ethylene glycol monostearate 8, 28

bis(2-Ethylhexyl) adipate 9

2-Ethylhexyl epoxy tallate 7, 29

2-Ethylhexyl sebacate 29

bis(2-Ethylhexyl) phthalate 8, 9

bis(2-Ethylhexyl) phthalate-3,4,5,6-d4 6

bis(2-Ethylhexyl) phthalate (DEHP) 2

bis(2-Ethylhexyl) terephthalate 4, 29

2,2'-Ethylidene-bis(4,6-di-tert-butylphenol) 15

N-Ethyl o,p-toluenesulfonamide 8, 29

F

F-300, F-1000, F-1500, F-2000, F-3000 22

Firemaster BP4A 24

G

Glycerol monooleate 8, 29

Glycerol monostearate 8, 29

Glycerol triacetate 7, 29

Glyceryl (triacetyl) ricinoleate 8, 29

H

Halowax 1013 25

Halowax 1051 25

Halowax 1099 25

2,2',3,4,4',5',6-Heptabromodiphenyl ether 23

Hercoflex® 900 29

2,2',4,4',5,5'-Hexabromodiphenyl ether 23

2,2',4,4',5,6'-Hexabromodiphenyl ether 23

Hexyl-2-ethylhexyl phthalate (Tech Mix) 7

Hi-Point PD-1 29

2-(2-Hydroxy-3,5-di-tert-amylphenyl) benzotriazole 15

I

Irganox® 245 16

Irganox® 259 16

Irganox® 565 16

Irganox 1035 40

Irganox® 1035 16

Irganox® 1081 16

Irganox® 1098 16

Irganox® 1425 WL 16

Irganox® 3114FF 16

Irganox® 3125 16

Irganox® E 201 17

Irganox® MD 1024 17

Isobutyl benzyl phthalate 7

Isobutyl cyclohexyl phthalate 7

Isonox® 132 17

Isonox® 232 17

Isopropyl isosteareate 8, 29

Isopropyl myristate 8, 29

Isopropyl palmitate 8, 29

J

Jayflex® 77 30

Jayflex® DIDP 30

Jayflex® DINP 30

Jayflex® DTDP 30

Jayflex® L11P-E 30

Jayflex® TINTM 30

K

Kemamide® E ultra 34

L

Laurex® 30

Lowinox® AH25 17

Lowinox® CPL 17

Lowinox® TBM-6 17

M

Markstat® 51 30

Markstat® 60 17

4,4'-Methylenebis(2-chloroaniline) 41

Methyl O-acetylricinoleate 30

Methyl oleate 8, 30

Mono-2-ethylhexyl phthalate 5

Mono-2-heptyl phthalate 5

Monobenzyl phthalate (mBzP) 4, 8

Monobutyl phthalate (mBP) 4, 8

Monocyclohexyl phthalate (mBP) 4

Monoethylhexyl phthalate (mEHP) 8

Monoethyl phthalate (mEP) 5, 8

Monoethyl phthalate 5

Monoisobutyl phthalate 5

Monoisononyl phthalate 5

Monoisopropyl phthalate 5

Monomethyl phthalate 5, 8

Mono-n-pentyl phthalate 5

Monoctyl phthalate 5

Morflex® 150 30

Morflex® 190 30

Morflex® 560 31

Morflex® x-1125 31

N

2-Naphthylamine 41

Naugard® 412S 17

Naugard® 445 17

Naugard® 635 17

Naugard® 956 18

Naugard® A 18

Naugard® B-25 18

Naugard® BHT 18

Naugard® HM-22 18

Naugard® J 18

Naugard® NBC 18

Naugard® PANA 18

Naugard® PHR 18

Naugard® PS-30 18

Naugard® PS-35 18

Naugard® Q Extra 18

Naugard® RM-51 19

Naugard® Super Q 19

Naugard® XL-1 19

N,N'-Dibutylthiourea 14

N,N'-Diethylthiourea 14

O

n-Octyl n-decyl phthalate (Tech Mix) 7

4,4'-Oxydianiline 41

P

Paraplex® G-30 31

2,2',4,4',5-Pentabromodiphenyl ether 23

2,2',4,4',6-Pentabromodiphenyl ether 23

Pentyl isopentyl phthalate 7

Perkacit® DPG 22

Perkacit® MBT 22

Perkacit® MBTS 22

Perkacit® NDBC 22

Perkacit® ZDEC 22

Plasthal® DINP plasticizer 31

Plasthal® ESO 31

Polycizer® butyl oleate 31

Polycizer® DP 500 31

Polyethylene glycol 200 dibenzoate 8, 31

Propylene glycol monostearate 8, 31

Propylene glycol ricinoleate 8, 31

Propyl gallate 19

n-Propyl oleate 8, 31

R

Resimene® 3520 22

Retarder AK 35

S

Santanox R 40

Santicizer® 141 31

Santicizer® 148 31

Santicizer® 160 32

Santicizer® 261 32

Santicizer® 278 19, 32

Santoflex® 6PPD 11

Santoflex 77PD 11

Santoflex IPPD 11

Saytex® 8010 25

SF100 11

Silquest® A-137 21

Silquest® A-187 21

Silquest® A-1100 21

Silquest® A-1102 21

Silquest® A-1289 21

Silquest® A-2171 21

Stearic Acid RG (rubber grade) 35

Stearic Acid TP 35

T

p-Terphenyl 25

m-Terphenyl 25

o-Terphenyl 25

2,2',4,4'-Tetrabromodiphenyl ether 23

Tetradecacloro-m-terphenyl 25

Tetradecacloro-o-terphenyl 25

Tetradecacloro-p-terphenyl 25

Tetrahydrofurfuryl oleate 8, 32

3,3',5,5'-Tetramethylbenzidine 41

bis(Thiobenzoyl) disulfide 35

4,4'-Thiodianiline 41

2,2'-(2,5-Thiophenediyi)bis(5-tert-butylbenzoxazole) 19

Tinuvin P 40

Tinuvin® PED 37

o-Tolidine 41

o,p-Toluenesulfonamide 8, 32

Tri-butoxyethyl phosphate 32

Tributyl phosphate 32

Tricapryl trimellitate 8, 32

Triethyl phosphate 32

Triisodecyl trimellitate 8, 32

2,2,4-Trimethyl-1,3-pentanediol-diisobutyrate 7

2,2,4-Trimethyl-1,3-pentanediol-isobutyrate 33

Trimellitate 33

2,4,4-Trimethylaniline 41

Tri(n-octyl, n-decyl) trimellitate 8, 32

tris(2-Butoxyethyl) phosphate 8

U

Ultranox 626 40

Ultranox® 626 19

Uvinul® 3000 37

Uvinul® 3008 37

Uvinul® 3040 37

Uvinul® 3049 37

V

Vinsol® powder 33

Vinsol® resin 33



Matériaux de Référence certifiés et Equipements de Laboratoire
www.techlab.fr Tél. (33) 03 87 75 54 29